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www.degiorgi.it e-mail info@degiorgi.it

INDICE

1 - SAFETY DIRECTIONS	4
2 - TECHNICAL FEATURES.....	6
3 - DEVICES OF CONTROL	6
4 - GENERAL INDICATION ABOUT THE UNIT	7
4.1 - CONNECTION OF THE EXHAUSTER	8
5 - INSTALLATION OF THE UNIT	9
5.1 - INSTALLATION	9
5.2 - ELECTRIC CONNECTION.....	9
6 - MOVEMENT	10
7 - GENERAL DIRECTIONS FOR USE	11
7.1 - INDICATION ABOUT THE USE OF THE UNIT	11
7.2 - WORKING OF THE MODEL DOWEL-PINS	12
7.3 - WORKING OF THE MODEL ACCU-TRAC.....	12
7.4 - WORKING OF THE MODEL TRICODENT	13
7.5 - SUGGESTED CUTTING TIME TO OBTAIN A CORRECT MOVEMENT OF THE TOOL	13
8 - PERIODICAL MAINTENANCE BY THE OPERATOR.....	13
8.1 - MAINS SUPPLY CORD	13
8.2 - ELECTRO-MAGNETIC PLANE	13
8.3 - LUBRIFICATION.....	14
8.4 - CLEANING OF THE FILTER.....	14
8.5 - CLEANING OF THE EXHAUSTER'S DEVICE.....	14
8.6 - CHANGING OF THE DISK.....	15
8.7 - CHANGING OF THE LAMP FOR ENLIGHTENED	15
9 - TROUBLESHOOTING GUIDE.....	16
10 - ACCESSORIES.....	17
11 - MAKER'S LIMITED GUARANTEE	18
12 - ELECTRIC DIAGRAMS.....	19
12.1 - ELECTRIC DIAGRAM.....	19
12.2 - ELECTRO-MAGNET POWER CIRCUIT	20

CE DECLARATION OF CONFORMITY
Directive 891392/CEE and 891336 and successive changes

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declares, under its own responsibility, that the apparatus:

Type: **ELECTRIC SAW**

described as:
electric saw making cut of plaster models

is in conformity with what is established by the directives
EN50081-1, EN55014, EN50082-1, EN292, 89/392/EEG,
modified 93/68/EEG and 91/368 EEG and 93/44EEG,
EMV89/336/EEG modified 92/31/EEG.

Signed: De Giorgi Ariberto
Managing Director

1 - SAFETY DIRECTIONS

WARNING

Look out with the use of electric tools and to avoid contacts with electric power, injuries and danger of fire; the following safety directions must always be observed.

Read and observe these indications before using the unit and keep with care these directions.

1. **Keep your working place in order** - Disorder on the working place may cause accidents.
2. **Enviromental conditions** - Don't expose electric tools to the rain. Don't use electric tools in wet or soaked places. Take care of good light conditions. Don't use electric tools near inflammable liquids or gas.
3. **Avoid electrification** - Avoid contacts with earthing objects, for ex.tubes, radiators, refrigerators.
4. **Keep children away** - Don't let third persons touch the unit or the cord, keep them far from the working place.
5. **Put your tool in a safe place** - Put your unit in a dry and sale place in order to not be accessible by children.
6. **Don' t overload your tool** - You'll work better and safer using the rating plate's power.
7. **Choice of the unit** - Don' t use low power tools or accessory that must be used with high power. Don' t use tools for work and purposes to which they aren't designed.
8. **Work attire** - During working, always dress properly, avoiding the use of loose-fitting dresses or too large sleeves.
9. **Use safety glasses** - Always use protective glasses. If you should work or something that produces dust, use a mask.
10. **Pay attention to the cord** - Don't lift the tool by the cord and don't use it to pull off the plug from the tap. Avoid contacts between the cord and high temperatures, oils, edges.
11. **Fix the piece fastened on which You are working** - Be sure that the piece is locked by a vice or in any case well blocked.
12. **Avoid unsafe positions** - Be sure to stay in a safe position, in order to maintain a good balance.

13. **Take care of your tools** - In order to work in a safe way, sharper and clean your tools, follow the assistance prescriptions and the directions about the change of using tools.
Often control the cord and, in case of damages, have qualified people change it. Control regularly the extension cords and change them if they are damaged; the side handgrips must be dry and without oil or fat.
14. **Changing of working tools** - Before making repairs and changing of tools, (f.e. saw-blades, points and in any case tools), always pull off the plug from the power tap.
15. **Don' t leave service keys on the unit** - Before using the tool be sure that service keys are taken away.
16. **Avoid unworked startings** - Before Your put the plug in the power tap be always sure that the switch is off.
17. **Don't use the unit in open places** - The unit must be used in well lighted and only indoors
18. **Always be careful working** - While working, don't be distracted. A lack of attention can cause accidents. In presence of other people in the working place, we recommend to maintain a suitable distance from the operator and the unit, in order to avoid unwanted and unintentional contact with machine parts in motion or residues of cuts.
19. **Always control that the tool isn' t damaged** - Before using the unit, check with attention its efficiency and the perfect working of the safety devices and if necessary of damaged parts.
Control the efficiency of the movable parts, that they aren' t blocked, that there aren' t broken parts, that the others are assembled in a right way and that all the conditions that could influence the correct efficiency of the unit are optimal. Safety devices or damaged parts must be repaired or changed by a Authorized Center of Assistance if different indications aren' t given.
20. **Whatch out!** - For your safety use only tools or accessories indicated in the direction for use or offered in proper catalogues. The use of different accessories or using tools, or anyway not indicate in the direction for use or in the catalogue, may cause accidents.
21. **Have qualified people change the tools** - This electric unit is in conformity with the safety directions in force. The repair of electric tools must be done only by authorized centers of assistance.
22. **Magnetic field** - The unit is endowed with a magnetic working plane which emits a significant magnetic field during the working phases. For this reason, people with medical devices (such as pacemakers) must keep a suitable distance from the unit.

2 - TECHNICAL FEATURES

		700/00
Voltage	V	220
Frequency	Hz	50
Base	mm	210/440
Height	mm	260
Weight	Kg	12
Blade	mm	80x16x0,2

3 - DEVICES OF CONTROL

1 **STARTING SWITCH** *Startin switch*

It's the switch which gives tension to the unit.

2-3 **Starting engine push-button**

Pushing both and simultaneously these push-buttons the disk begins its rotation. While working they must always be, because the release of only one causes the disk to stop immediately. They also control the start of the exhauster and the activatier of the electro-magnet.

4 **Fuseholder circuit 220V**

It contains the fuse that protects the chief circuit at 220V.

5 **Fuseholder circuit 220V**

It contains the fuse that protects the exhauster's engine or the one of the external exhauster.

6 **Fuseholder circuit 24V**

It contains the fuse that protects the low tension circuit.

7 **Screw for the locking model**

The locking of this screw blocks the plaster model on the modelholder.

8 **Lever for the locking of the inclination**

The achivation on of this lever allows the locking of the modelholder in the position preventively chosen by the operator.

9 **Warning light**

The putting on of this light is got pressing the two push-buttons of the starting engine; it points out that the electro-magnet is on. This light will always be lighted while working, it will be put out only at the end of working and a few seconds later than the engine. The electro-magnet will be on for the necessary time to avoid mistakes and danger caused by faulty operation when the disk is still in movement of inertia.

4 - GENERAL INDICATION ABOUT THE UNIT

The saws are units built expressly to make the cut of plaster models in order to separate the dies. They are formed of a housing made of a painted steel plate on which are fixed the following main assemblies.

☐ MOVING ARM

Includes the square section that acts as structure, the engine, the members of transmission and the tool for working.

☐ ELECTRO-MAGNETIC PLANE

Formed by the engine, the filter, the tube of connection and the dusts' conveyor.

☐ ELECTRIC SYSTEM

Includes the switches, the push-buttons and the printed circuit.

In order to increase the safety of these units there are the following protections:

- a bladecover that avoids the casual contact with the disk in the upper zone;
- a protection made of plastic-transparent material with the function of holding bodies or particles that could project themselves towards the operator during the working.

TO THIS PURPOSE THE OPERATOR MUST ALWAYS WEAR PROTECTIVE GLASSES AND A DUST MASK;

- an electro-magnetic plane endowed with a high power that guarantees a great stability to the modelholder base.

It generates a magnetic field during its operation, so don't bring near to it materials that could be damaged by it, such as magnetic tapes, floppy disks, etc. In particular, people with medical devices, such as pacemakers, must be careful;

- a returning spring of the moving arm that brings it back to its stopping position;
- a driving cog belt that, besides being free from maintenance, allows the tool to stop itself in a short time after the release of the push-buttons;
- a two-hands command that necessitates suppression of both the two push-buttons in order to activate and to maintain the disk in rotation. They are crossed by a low electric tension.

As mentioned earlier, these units are projected to cut plaster models, for this reason it's important:

- to cut always dry plaster, because wet plaster may pug itself with the breaking of the disk;

- don't use materials different from the ones allowed;
- connect an exhauster to the electric saw where such is lacking using the connections placed at the back of the unit (See par. 4.1);
- always in good conditions the system of suction cleaning often the filter (see the chapter about the maintenance in this handbook par. 8.4).

4.1 - CONNECTION OF THE EXHAUSTER

This operation must be done by qualified personnel.

Important thing: before connecting the exhauster with the unit, be sure that it hasn't a capacity lower than 15 lt/sec and that the highest power isn't greater than 800 W.

This condition must be observed in order not to hinder the working of the unit.

In order to make the connection, follow these indications:

- 1 pull off the plug from the power tap of the unit;
- 2 connect, through a flexible tube the exhauster to the unit plugging it in its proper tube. The tube must have an inner diameter longer than 20 mm. in order to guarantee the passage of air at least equal to 3 cmq., and a length (according to your necessity) so short in order to reduce the loss of the load. In addition it mustn't have excessive crushings or foldings;
- 3 connect the main supply cord of the exhauster to the power tap on the back side of the unit.

If the plug is different from the required model, change it with the one that is supplied with the unit.

After the connection, the exhauster will start operating simultaneously with the rotation of the disk, putting on the push-buttons.

5 - INSTALLATION OF THE UNIT

5.1 - INSTALLATION

- Remove the unit from the package by at the opening side taking care of not to collide against the disk.
- Check that all the parts of the unit are whole. If this should not be the case, immediately contact the retailer.
- Lay the unit on a fixed plane in a well lighted and dry place and, free with space of at its sides, at least 30 cm.
- Now the unit is assembled, ready to work.
- Connect the unit to the exhauster if such should be lacking.
- Don't throw the package away, but conserve it for next movements.

5.2 - ELECTRIC CONNECTION

- Before connecting the unit be sure that the tension of service indicated on the data plate agrees with the one of the system and that it is endowed with earth wire.
- Connect the plug of the electric flying cord coming from the unit to the proper tap. Make this operation when the general switch is off.
- The main supply of the unit must be done according with the instructions of CARLO DE GIORGI srl that, however, isn't responsible the connection fax.

The endowed mains supply cord is of the model H05 with a section consistent with the highest absorbed power (indicated in the technical features).

The supplied cord has a lenght of 2 m. In case of using extensions, be sure that:

- the collocation of the cord doesn't cause damages ;
- the point of connection to the eletric power font garanteos the presence of all the necessary protections for the directions in force.

6 - MOVEMENT

Conserve the package of the unit that will be useful if you have to move it to a different working place.

In this case, follow the operations indicated below. We suggest you dismantle the disk in order to avoid it being accidentally damaged casual damage.

- A) Switch off the unit, disconnect the electric plug and later pull it out of the unit.
- B) Put the modelholder base in its proper box with the supplied keys, and put all in the pack in its place.
- C) Put the unit in its box in its proper place.
- D) Put in the electric plug and the handbook of the directions for use and maintenance, then put on the cover and close the box.

7 - GENERAL DIRECTIONS FOR USE

7.1 - INDICATION ABOUT THE USE OF THE UNIT

The operator must be in front of the unit in order to use his right hand on the moving arm and on its push-button, while with his left hand he must activate the controls placed on the side of the base.

Concerning the use of the unit, follow these operations chronologically:

- 1 fix the model on the modelholder being sure that it can't be moved during working;
- 2 lay the complex base+model on the electro-magnet plane;
- 3 swith on the unit using the push-button, simultaneously also the light of working will be switched on;
- 4 put the complex base+model on the working plane getting the disk to cut in the place desired.

For this purpose you can change the inclination of the modelholder using the proper lever;

WARNING

Before start to work, be sure that the protection screen is led down.

- 5 push the push-buttons that make the disk move and the exhauster simultaneously, they also command the operating of the magnetic piano that fixes the base on it. The warning light will draw attention to this situation;

WARNING

Never set in motion the unit when the disk is on the model: when it starts its rotation it must be free.

- 6 go on pushing the starting push-buttons let down the moving arm and make the cut on the model until its achieved the wanted depth.

WARNING

More pressure doesn't mean necessarily a quick cut, on the contrary it can cause an excessive heating of disk and unit.

- 7 at the end of the cut put the disk back in its relaxed position and pull away from the starting push-buttons;
- 8 if you will not use the unit within a short time, switch off the unit and pull out the plug from the power tap.

With the use of the electric saws it is possible to work the plaster model by three different systems that is to say Dowel-Pins (double casting with supports), Accu-trac and Tricodent using each time the modelholder made specially.

7.2 - WORKING OF THE MODEL DOWEL-PINS

This modelholder is supplied with the unit. In order to work the model, it is very important to insert the supports in the impression in a vertical position (in this case would be excellent the use of a pinning unit).

Therefore be careful not to make cuts that can interfere with the support or part of it, because this situation involves the quick wear of the disk or even the break of some of its sectors, in this case exchange it immediately with a new one.

- The right position of the model is the one below the burholder shaft.
- It's also important, during the cut that the disk doesn't interfere with the screws on the modelholder, these must be regulated on the lower position compatible with a good fixing of the model.

7.3 - WORKING OF THE MODEL ACCU-TRAC

The modelholder ACCU-TRAC is an accessory on request. This kind of model must be put on the modelholder in a way that the teeth of the mechanism, securing the taking catch themselves perfectly in the knurled low part of the model.

In this way it will be able to be blocked, pressing the knob in any positions at the condition that the three mechanisms are always in taking position.

Thanks to this system, the zone of the disk concerning the working is the one in front of the burholder shaft this allows us to make use of a greater depth of cut.

- Place than the complex model+base shifted forward as regards the disk.
- Be careful that the disk doesn't hurt the modelholder or its mechanisms during the working.
- The plate of support of the base must never come out from the central core of the electro-magnet because in this case the traction will be nothing.
- In order to avoid the fall of the sectioned dies, pouring also the inner zone of the model that will take them in their position.

7.4 - WORKING OF THE MODEL TRICODENT

The modelholder TRICODENT is an accessory on request.

It is formed by a sort of terminal that blocks the model on its internal part: it's consequently necessary pouring it at all.

Also this kind of model agrees with the before sayd speech according to which the zone of the disk concerning the working is the one in front of the burholder shaft.

- Put the complex base+model shifted forward as regard to the disk.
- Be careful that the disk doesn' t hurt the modelholder.
- The plate of support of the base must never come out of the central core of the electro-magnet because in this case the traction will be nothing.

7.5 - SUGGESTED CUTTING TIME TO OBTAIN A CORRECT MOVEMENT OF THE TOOL.

To make the cut on a pouring model with the system Dowel-Pins according to the height and the section of the model, the working time can change from a minimum of 4 seconds to over 10 seconds.

Instead the models ACCU-TRAC and TRICODENT have longer working time, because the section of the model includes also the base, so the cutting time increases from a minimum of 10 seconds to a maximum of 18-20 seconds.

It's important then to use a slow movement that allows to empty the asported plaster, obtaining a preciser and straighter cut.

8 - PERIODICAL MAINTENANCE BY THE OPERATOR

All the operations for maintenance must be done when the unit is stopped, disconnected by the mains supply.

8.1 - MAINS SUPPLY CORD

Control periodically the conditions of the power cord. If it is damaged, make qualified people change it.

8.2 - ELECTRO-MAGNETIC PLANE

Take always the electro-magnetic plane clean, if necessary, grease it with vaseline oil to protect it from oxidation.

8.3 - LUBRIFICATION

The unit is endowed with waterproof bearings, so it isn't necessary the lubrication of the moving parts.

8.4 - CLEANING OF THE FILTER

Disassemble the sponge filter each 5 or 6 days in which the saw was used.

In order to do that, unscrew the locking knob, than take out the filter and clean it with compressed air, later put it such a before. The obstructed filter causes a wrong aspiration of the dusts produced during the working. If the filter was damaged, change it asking only to our technical assistance.

8.5 - CLEANING OF THE EXHAUSTER'S DEVICE

After a long working time it can be necessary eliminate the plaster sediments that are formed inside the exhauster's devices, for this reasons their section can tighten so much that compromises the right working of the exhauster.

In order to make this operation it's necessary to disassemble the side where the dusts are conveyed after taking off the screws that fix it.

8.6 - CHANGING OF THE DISK

The disk must be changed each time its conditions doesn't allow of using it in absolute safety conditions.

Change it without delay if it appears bent or if one or more sectors are lacking or bent in consequence of impacts. The wear of the disk isn't easily identified only with a visual check, anyway its deterioration is identified valuing the operator's effort during the working. In fact a wearing disk lacks of its abrasive power.

Don't assemble absolutely disks made of HSS (super-quick steel), because they chip and break shattes all of a sudden.

WARNING

Before changing the disk switch off the unit and pull off the plug from the power tap.

Put the keys of 4 mm in the hexagon of the sinistrorse screw that blocks the disk, then insert the support in the hole made on the shaft that in this way will be blocked.

Unscrew the screw in a clockwise direction as far as the disk will be free. Reassemble the new disk and make the before said operations in the contrary direction.

The disk hasn't an its own direction of rotation, so it can be assembled in both the directions.

8.7 - CHANGING OF THE LAMP FOR ENLIGHTENED

In the case you must change this lamp, use another with the same features. In order to take it away from the lampholder it's enough a little traction towards the inside and the contacts will be open making it free.

Instead to reassemble it you must lay and press it in the direction of the lampholder.

9 - TROUBLESHOOTING GUIDE

All this operations must be done when the unit is disconnected. Don't make operations of maintenance regulations or repairs that aren't indicated in this handbook for use and maintenance. In case of necessity contact our technical assistance.

TROUBLES	CAUSES	SUGGESTIONS
The unit doesn't work.	The electric power doesn't arrive to the unit.	Be sure that the cord of connection between the unit and the electric system is correctly connected in the respective taps. Control the respective fuse (for this operation turn to qualified personal).
The lamp for enlightened doesn't work.	Break of the filament. Control the respective fuse that must be burn.	Take it off from the lampholder and be sure that it's integral, on the contrary case change it with a new one. Change the respective fuse (for this operation turn to qualified personal).
The plaster of the model has a tendency to clog and to stop the disk.	The model could have a too high level of dampness.	Control the drying process of the model.
During the working the disk has a tendency to vibrate or to break the model or anyway doesn't work correctly.	The disk is damaged. The disk is assembled in a wrong way or extraneous particles could be between it and the planes for the locking. The disk is worn.	Control that it hasn't had impact, that it hasn't bent or lacking parts. In this case change it. Disassemble it and clean its surfaces. (see par. about MAINTENANCE) . Change the disk (see par. about MAINTENANCE).

If the imperfection isn't resolved although the made checks, we advise you to ask the retailer or the maker or all the damages connected with: transmission, engine or electric parts you must ask the maker.

10 - ACCESSORIES

ENCLOSED ACCESSORIES:

- Electric cord with a tap EN 60320-C13 and a plug SCHUKO CEE7 -S.S.VII.
- Setscrew wrench of 4 mm.
- Cylindrical plug with a diameter of 3 mm.
- Rotary base.
- Modelholder DOWEL-PINS.
- Lamp and fuse of replacement kit.
- Flying plug.

ACCESSORIES DISPOSABLE ON REQUEST

- Modelholder ACCU-TRAC.
- Modelholder TRICODENT.
- Diamond disk 80x0.25 mm. (hole 16 mm).

11 - MAKER'S LIMITED GUARANTEE

The guarantee is limited CARLO DE GIORGI srl guarantee that the saw operated by electric engine is free from imperfections of material and processing for a year from the purchasing date.

The guarantee is valid only if the product was bought from CARLO DE GIORGI srl or authorized retailers.

This guarantee isn't valid if the product was used in an illicit way, if was altered, changed or repaired by unauthorized personal by CARLO DE GIORGI srl or by unauthorized services' agents.

In the case that the product shows, in one year from the purchasing, imperfections of material or processing it must be sent to CARLO DE GIORGI srl or to an authorized by CARLO DE GIORGI srl service's centre, paying the travelling expenses, enclosing name, address, the purchasing date's proof and a short declaration of the noticed imperfection.

CARLO DE GIORGI srl as it likes, well make the repair or the change of the imperfect product without any expense by the purchaser.


Possible repairs or changes are guarantee in the above described way for the rest of the original guarantee time. The guarantee is limited to the repair and change of the imperfect product.

(There are no other implicit or explicit guarantees and in addition CARLO DE GIORGI srl will be not responsible of casual or special damages or other damages, costs or expenses except the ones of repair or change as above described).

12 – ELECTRIC DIAGRAMS

12.1 - Electric diagram

carlo de giorgi s.r.l.

Via Tonale, 1 – 20021 Baranzate di Bollate (MI) – Italy
Tel. +39.02.356.15.43 r.a. – Fax +39.02.356.18.08
<http://www.degiorgi.it> E-Mail: info@degiorgi.it