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SUMMARY

1 - INTRODUCTION.....	4
1.1 - FOREWORD	4
1.2 - DESCRIPTION OF THE DEVICE	5
1.3 - INTENDED USE	5
1.4 - SAFETY REQUIREMENTS.....	5
2 - IDENTIFICATION DATA.....	8
2.1 - IDENTIFICATION DATA	8
2.2 - DATA PLATE OF THE DEVICE.....	8
2.3 - DATA PLATE OF THE SCALER HANDPIECE	8
3 - TESTING OF THE DEVICE.....	9
4 - TRANSPORTATION	9
5 - LIST OF CONTENTS	9
6 - INSTALLATION	11
6.1 - SAFETY REQUIREMENTS DURING INSTALLATION	11
6.2 - INITIAL SETUP.....	12
6.3 - CONNECTION OF THE ACCESSORIES.....	12
7 - CONTROLS	14
7.1 - DESCRIPTION OF THE CONTROLS	14
7.2 - DESCRIPTION OF THE DISPLAY AND FUNCTIONS	16
7.3 - SAFETY REQUIREMENTS DURING USE.....	16
7.4 - PROTECTION SYSTEMS AND ALARMS	18
7.5 - INSTRUCTION FOR USE	19
7.6 - RULES FOR KEEPING THE DEVICE IN PROPER WORKING ORDER	19
7.7 - SETTINGS PERMITTED ACCORDING TO INSERT TYPE	20
8 - CLEANING, DISINFECTION AND STERILIZATION	20
8.1 - CLEAN FUNCTION - CLEANING OF THE LIQUID CIRCUIT	20
8.2 - CLEANING AND DISINFECTING THE CASING OF THE APPARATUS	21
8.3 - STERILIZATION PROCEDURE	22
8.4 - AUTOCLAVE STERILIZATION OF THE HANDPIECE	22
8.5 - AUTOCLAVE STERILIZATION OF THE INSERTS	24
8.6 - AUTOCLAVE STERILIZATION OF THE WRENCH FOR TIGHTENING THE INSERTS	24
8.7 - AUTOCLAVE STERILIZATION OF THE PERISTALTIC PUMP TUBE	24
8.8 - AUTOCLAVE STERILIZATION OF THE CONNECTION BETWEEN CORD AND PERISTALTIC PUMP TUBE CONNECTION	25
8.9 - AUTOCLAVE STERILIZATION OF THE HANDPIECE SUPPORT	25
9 - REGULAR MAINTENANCE.....	25

10 - REPLACEMENT OF THE FUSES..... 26

11 - DISPOSAL PROCEDURES AND PRECAUTIONS 26

12 - TIPS 27

12.1 - SHARP INSERTS..... 27

12.2 - SMOOTHING INSERTS 27

12.3 - BLUNT INSERTS..... 27

13 - SYMBOLS 28

14 - TROUBLESHOOTING..... 29

15 - TECHNICAL DATA..... 32

16 - AFTER SERVICE..... 33

17 - ENVIRONMENTAL PROTECTION 33

18 - MANUFACTURER’S RIGHT 33

19 - PLEASE CONTACT, FOR TECHNICAL DATA 33

20 - DECLARATION OF CONFORMITY..... 33

20.1 - PRODUCT CONFORMITY THE FOLLOWING STANDARDS 33

20.2 - DECLARATION OF CONFORMITY - EMC..... 34

21 - GUARANTEE 38

1 - INTRODUCTION

1.1 - FOREWORD

Before proceeding with the installation, use, maintenance or any other activities on the equipment please read the manual carefully.

Important: To avoid causing personal injuries or damages to property, read all the points concerning “safety requirements” contained in this manual with particular attention.

Depending on the level of risk involved, safety requirements are classed under the following indications:

Danger: (always referred to personal injury).

Warning: (referred to possible damage to property).

The purpose of this manual is to ensure that operators are aware of the safety requirements, of the installation procedures and of the instructions for the correct use and maintenance of the device.

The user is not authorized to tamper with the equipment under any circumstances.

If any problems are encountered, please contact a Carlo De Giorgi Service Centre.

Any attempts on the part of the user or any unauthorized personnel to tamper with or alter the device will invalidate the warranty and release the Manufacturer from any liability in respect of any harm or damage to persons or property.

The information and illustrations contained in this manual are updated to the date of publication indicated on the last page.

Carlo De Giorgi is committed to continuous updating of their products, which may entail some changes to components of the equipment. If there are any discrepancies between the descriptions contained in this manual and your equipment, please contact your dealer or the Carlo De Giorgi after-sale service for explanations.

Using this manual for purposes other than those relating to the installation, use and maintenance of the equipment is strictly prohibited.

1.2 - DESCRIPTION OF THE DEVICE

Thanks to its controlled three-dimensional ultrasound oscillations, the original Ultrasurgery technique brings in a new age for osteotomy and osteoplasty in Implantology, Periodontology, Endodontics and Orthodontic Surgery. Its main features are:

- Micrometric cutting: Maximum surgical precision and intra-operative sensibility;
- Selective cutting: Maximum safety for the soft tissues;
- Cavitation effect: Maximum intra-operative visibility (bloodless field);

The equipment has an automatic tuning circuit that offsets the wear of the inserts, thus ensuring work in constant conditions of maximum efficiency.

1.3 - INTENDED USE

The Ultrasurgery is a piezoelectric device for bone surgery that allows osteotomy and osteoplasty techniques to be applied in almost any anatomical situation. This equipment can be used in the following fields: (1) oral surgery; (2) orthopaedic surgery; (3) maxillofacial surgery; (4) cosmetic surgery; (5) neurosurgery; (6) otolaryngology.

This equipment cannot function in places where there is an inflammable atmosphere (anaesthetic mixture, oxygen, etc.).

1.4 - SAFETY REQUIREMENTS

Carlo De Giorgi will not accept any liability for direct or incidental personal injury or damage to property in the following cases:

- (1) If the equipment is used for purposes other than that for which it is intended.
- (2) If the equipment is not used in accordance with all the instructions and requirements described in this manual.
- (3) If the wiring system in the room where the equipment is used does not comply with the application standard and appropriate requirements.
- (4) If any assembly operations, extensions, settings, alterations or repairs have been carried out by personnel not authorized by Carlo De Giorgi.
- (5) If the environmental conditions in which the device is kept and stored do not comply with the requirements indicated in the chapter on technical specifications.

Danger: Qualified and specialized personnel.

This equipment may be used only by specialized and suitably trained personnel such as surgeons. If correctly used, this equipment does not give rise to side effects. Improper use, on the other hand, will give rise to transmission of heat to the tissues.

Danger: Intended use.

Use the equipment solely for the purpose for which it is intended (see point 1.3). Failure to comply with this requirement could lead to serious harm to the patient and/or to the operator and/or damage to/failure of the equipment.

Danger: Contraindications.

Do not use the Ultrasurgery on patients with pacemakers or other implantable electronic devices. The same requirement applies also to the operator.

Danger: Contraindications.

An electrosurgical knife could interfere with correct functioning of the device.

Danger: Cleaning, disinfection and sterilization of new or repaired products.

All new or repaired products are delivered in no sterile conditions. Before being used for treatments, all new or repaired products should be cleaned, disinfected and sterilized following strictly the instructions provided under point 8.

Danger: Use only original Carlo De Giorgi accessories and spare parts.

Danger: Check the condition of the device before treatment. Always make sure that there is no water under the device. Before each treatment always check that the equipment is in proper working order and that the accessories are undamaged. Do not carry out the treatment if any problems are encountered while operating the device. If the problems concern the equipment contact an authorized technical service centre.

Danger: Breakage and wear of the inserts. The high-frequency vibrations may very occasionally lead to breakage of heavily worn inserts. Inserts of which the shape has been changed or which are otherwise damaged are liable to break during use. Any such inserts should absolutely not be used. It is necessary to instruct the patient to breathe through his nose during the treatment in order to avoid ingestion of the broken off fragment of the insert.

Danger: Do not install this equipment anywhere there is a risk of explosions. This equipment cannot function in places where there is a flammable atmosphere (anaesthetic mixture, oxygen, etc.).

Danger: Personnel injury. The footswitch of the Ultrasurgery must not be activated when the door of the peristaltic pump open. (Fig. 4 - Ref. B). Moving parts could injure the operator.

Danger: Contraindication. Do not carry out this treatment on metal or ceramic prosthetic artifacts. The ultrasonic vibrations could lead to serious damage to such artifacts.

Danger: Contraindication. After autoclave sterilization of the handpiece, wait for it to cool down completely before using it.

2 - IDENTIFICATION DATA

2.1 - IDENTIFICATION DATA

An exact description of the model including the serial number of the equipment will make it easier for our After-Sale Service to respond quickly and efficiently to your enquiry. Always provide the above information whenever you contact a Carlo De Giorgi Service Centre.

2.2 - DATA PLATE OF THE DEVICE

Each device has its own data plate (Fig. 1), on which technical specifications and serial number are indicated. The data plate is on the rear of the device. The remaining data are included in this manual (see point 15).


Fig. 1

2.3 - DATA PLATE OF THE SCALER HANDPIECE

The serial number of the Ultrasurgery handpiece is engraved on the metallic sheath (Fig. 2).


Fig. 2

3 - TESTING OF THE DEVICE

1. All the devices are checked and tested by Carlo De Giorgi completely, including all the parts.
2. When testing, all the parts will work in intermittent operation.
3. The test emphasizes all problems deriving from any faulty parts.
4. This procedure ensures the function and reliability of all the parts.

4 - TRANSPORTATION

1. Avoid any excessive shock or shaking, and keep covered during transportation.
2. Do not transport together with dangerous materials.
3. Avoid exposure to sunlight, rain, or snow during transportation.

5 - LIST OF CONTENTS

The content may vary in case of promotional campaigns.

Warning: Handpiece and cord cannot be detached.

Name	Quality	Ref
Peristaltic pump	1	Fig. 3 - Ref. A
Device	1	Fig. 3 - Ref. B
Torque wrench	1	Fig. 3 - Ref. C
Ultrasurgery handpiece complete with cord	1	Fig. 3 - Ref. D
Connection for the cord and tube of the peristaltic pump	1	Fig. 3 - Ref. E
Connector of handpiece	1	Fig. 3 - Ref. F
Plug of Footswitch	1	Fig. 3 - Ref. G
Footswitch	1	Fig. 3 - Ref. H
Insert kit	Marked on the packing list	Fig. 3 - Ref. I
Insert	Marked on the packing list	Fig. 3 - Ref. J
Input of power-supply cable	1	Fig. 3 - Ref. K
Output of power-supply cable	1	Fig. 3 - Ref. L
Surgical tray	1	Fig. 3 - Ref. M
Rod for supporting the bag	1	Fig. 3 - Ref. N
Support for the handpiece	1	Fig. 3 - Ref. O
Pump tube	1	Fig. 3 - Ref. P


Fig. 3

6 - INSTALLATION

6.1 - SAFETY REQUIREMENTS DURING INSTALLATION

Danger: The wiring system of the premises where the apparatus is installed and used must comply with the applicable standards and the relevant electrical safety requirements.

Danger: Do not install the apparatus in places where there is a risk of explosion. The apparatus may not be used in areas where there are flammable atmospheres (anaesthetic mixtures, oxygen, etc.).

Danger: Install the apparatus in a place where it will be protected from blows and from accidental sprays of water or other liquids.

Danger: Do not install the device on or near sources of heat. Install it such a way that there is an adequate circulation of air around it. Leave sufficient free space around it, in particular with reference to the fan on the rear (Fig. 5).

***Danger: Personnel injury.
The footswitch of the Ultrasurgery must not be activated when the door of the peristaltic pump open. (Fig. 4 - Ref. B).
Moving parts could injure the operator.***

Warning: Do not expose the device to direct sunlight or to sources of UV light.

Warning: The device is transportable; however, it must be handled with care when it is moved.

Warning: Before connecting the power supply cord to the device, make sure that the electrical contacts are perfectly dry. If necessary, dry them with an air syringe.

Warning: To avoid risks of electric shock, only connect this device to a supply line equipped with earthing systems.

6.2 - INITIAL SETUP

The device must be installed in a suitable environment for its use:

- install the device following the instructions, starting from chapter 6 "Installation", referencing the images contained in this manual;
- photocopy the purchase invoice, showing date and serial number of the device to grant warranty and traceability of equipment itself;
- keep a copy of the invoice with the device, and send it together with the device back to technical assistance, in case of repairs.

6.3 - CONNECTION OF THE ACCESSORIES

The accessories listed should be connected with the Ultrasurgery following the process here described.

- 1 Insert the silicone tube into the peristaltic pump, proceeding as follows:
 - Open the lid (Fig. 4 - Ref. A) as far as it will go;
 - Position the tube on the impeller (Fig. 4 - Ref. B).
 - Close the lid completely (Fig. 4 - Ref. C).


Fig. 4


Danger: Personnel injury.

The footswitch of the Ultrasurgery must not be activated when the lid of the peristaltic pump open. (Fig. 4 - Ref. B). Moving parts could injure the operator.

- 2 Insert the rod for supporting the bag into the holes provided for it (Fig. 5 - Ref. A);
- 3 Connect the footswitch to the main body of the device by inserting the plug into the footswitch socket (Fig. 5 - Ref. D);
- 4 Plug the power cable into the connector on the main body of the device (Fig. 5 - Ref. E) and then into the power outlet;

- 5 Insert the handpiece support into the two holes provided for it, then plug the handpiece cord in the main body connector according to the red dots (Fig. 5 -Ref. C; Fig. 6);
- 6 Insert the handpiece tube in the pump connector (Fig. 3 - Ref. P);
- 7 Put the handpiece on the support (Fig. 3 - Ref. O);
- 8 Connect end of the tube of the peristaltic pump;
- 9 Connect the flow-control system to the bag containing the appropriate liquid for the treatment.
- 10 Use the torque wrench to screw the insert (Fig. 7) till the clattering noise.
- 11 press the “on/off” button (Fig. 5 - Ref. B); the device can be used.

Fig. 5**Fig. 6****Fig. 7**

7 - CONTROLS

7.1 - DESCRIPTION OF THE CONTROLS

This section illustrates the parts of the front panel of the Ultrasurgery unit, enabling the controls described in this manual to be located immediately.

(1) *Description of bone function:*


Fig. 8

(2) Description of root function:


Fig. 9

(3) Description of clean function:


Fig. 10

7.2 - DESCRIPTION OF THE DISPLAY AND FUNCTIONS

There are three functions of bone, root, clean for this Ultrasurgery.

a) *BONE function (Fig. 8)*

In bone function, the power setting is unavailable, but it is possible to adjust the water setting. There are four bone settings, as follows:

- 1 Quality1: very high bone density.
- 2 Quality2: high bone density.
- 3 Quality3: middle bone density.
- 4 Special: low bone density.

b) *ROOT function (Fig. 9)*

In this function, both the water and power settings are available; two modes are available as follow:

- 1 Perio.
- 2 Endo.

c) *Clean function (Fig. 10)*

In this function, press the footswitch and hold it on for 3 seconds; the device will clean the tube from water automatically in 25 seconds.

7.3 - SAFETY REQUIREMENTS DURING USE

Danger: Contraindications.

Do not use the Ultrasurgery on patients with pacemakers or other implantable electronic devices. This requirement also applies to the operator.

Danger: Breakage and wear of the inserts.

The high-frequency vibrations and wear may, very occasionally, lead to breakage of heavily worn inserts. Inserts of which the shape has been changed or which are otherwise damaged are liable to break during use. Any such insert should definitely not be used. It is necessary to instruct the patient to breathe through his nose during the treatment in order to avoid ingestion of the broken fragment of the insert.

***Danger: Control of infections.***

For maximum safety of both the patient and the operator, clean, disinfect and sterilize the piezo electronic handpiece, the inserts and the torque wrench after each treatment.

***Warning: Contraindication***

Do not carry out this treatment on metal or ceramic prosthetic artifacts. The ultrasonic vibrations could lead to serious damage to such artifacts.

***Warning: Contraindication***

After autoclave sterilization of the handpiece, wait for it to cool down completely before using it.

Warning: The electrical contacts inside the cord connector must be dry.

Before connecting the handpiece to the device, make sure the electrical contacts of the connector are perfectly dry, in particular after the autoclave sterilization cycle. If necessary, dry the contacts by blowing air onto them with the air syringe.

Warning: To use the device correctly, it is necessary to press the footswitch and start it up without letting the insert rest on the part to be treated. This will allow the electronic circuit to detect the point where resonance of the insert it best without any interference, thus enabling optimum performance. If this is not done, contact with the part to be treated or with other surfaces before start-up could cause the activation of the protection systems.

Warning: for spray treatment, use only inserts that allow passage of liquids.

7.4 - PROTECTION SYSTEMS AND ALARMS

The device has a diagnostics circuit that is used to recognize activation of the protection system and of the alarms. These are shown on the display, as follows:

Warning code	Possible Cause
Warn 01	There is a signal transmission problem.
	Switch the device off and then on again, If the problem persists contact Carlo De Giorgi service centre.
Warn 02	Tuning circuit not working properly.
	Handpiece failure.
Warn 03	Fan failure.
Warn 04	Pump failure.
Warn 05	Voltage too high.
Warn 06	Insert not correctly secured to the handpiece or insert worn, broken or deformed.
	Handpiece failure.
Warn 07	Tuning circuit not working properly.
	Please restart the device, if the problem persists, stop using it and call the Carlo De Giorgi service centre immediately.

7.5 - INSTRUCTION FOR USE

- a) Open the air intake on the drip system.
- b) Screw the chosen insert onto the Ultrasurgery handpiece until it is flush against it.
- c) To use the torque wrench correctly (Fig. 7) proceed as follow:
 - 1 Hold the body of the handpiece firmly;

Warning: Do not grip the end part of the handpiece or the cord, only the plastic casting (Fig. 7) and do not turn it while fastening the insert in place;

- 2 Turn the wrench in a clockwise direction until the clutch engages (making a clicking sound);
 - 3 The insert is now properly tightened in place.
- d) Make sure that the Ultrasurgery handpiece is correctly connected to the handpiece connector (Fig. 5 - Ref. C, Fig. 6).
 - e) Check the display to see the type of power that has been set. If the type of power required is different from the type that has been set, use the "M" key (see 7.1) to switch to the correct one.
 - f) Check the display to see the power level that has been set, if the type of power required differs from the level that has been set, use the "+", "-" keys (Fig. 9) for selecting, depending on the type of function that has been set.
 - g) Check the display to see the delivery rate of the peristaltic pump, if the delivery rate required is other than the level that has been set, use the "+", "-" keys (see 7.1) to choose, depending on the type of function that has been set.

7.6 - RULES FOR KEEPING THE DEVICE IN PROPER WORKING ORDER

- 1 Check the state of wear of the inserts periodically and replace any for which a drop in performance is noted;
- 2 Do not alter the shape of the inserts by bending or filing them;
- 3 Replace any insert that has become deformed or damaged by impacts;
- 4 Always make sure that any threaded parts and their contact surfaces are perfectly clean;
- 5 If an insert becomes too worn, the device will stop working.

7.7 – SETTINGS PERMITTED ACCORDING TO INSERT TYPE

The following table shows the Function and Power settings permitted for each insert, for the correct use of the device.

Insert	Mode	Power
US1-US2-US3-US4-US5-US6-US1L-US1R	BONE	Quality1-Quality3, Special
UL1-UL2-UL3-UL4-UL5	BONE	Quality1-Quality3, Special
UC1	BONE	Quality1-Quality3, Special
UI1-UI2-UI7-UI8-UI9	BONE	Quality1-Quality3, Special
UP1-UP2-UP3-UP4-UP5-UP6-UP7	ROOT	Perio
UE1-UE2-UE3-UE4	ROOT	Endo

8 - CLEANING, DISINFECTION AND STERILIZATION

8.1 - CLEAN FUNCTION - CLEANING OF THE LIQUID LINE

Warning: failure to carry out cleaning of the tubes will lead to crystallization of salts that can seriously damage the equipment.

Warning: handpiece and cord cannot be detached.

- 1** Change the bag containing water (demineralized water is recommended);
- 2** Check whether the water system is connected correctly;
- 3** Start the **CLEAN** function (Fig. 10);
- 4** Press the footswitch to start the cleaning cycle, as soon as the peristaltic pump starts up, a status bar will appear on the display to indicate progressively the time remaining to completion of the **CLEAN** cycle. The cycle lasts for 25 seconds and cannot be stopped.

- 5 *Once the cleaning cycle has been completed, the device exits from the CLEAN function and returns to BONE function (Fig. 8);*
- 6 *On completion of the cleaning operations, empty the tubes and dry the accessories that have been through the cleaning cycle.*

8.2 - CLEANING AND DISINFECTING THE MAIN UNIT

***Danger:** The main unit is not protected against the penetration of liquids. Do not spray liquids directly onto the surface of the main unit.*

***Danger:** The device cannot be sterilized.*

After each treatment carry out the following operations:

- 1 *Remove the insert from the handpiece;*
- 2 *Clean and disinfect the surfaces of the main unit, the cords and their connectors using a cloth moistened with a mild detergent or disinfectant solution with a neutral pH (pH 7). Follow carefully the instructions given by the manufacturer of the disinfectant solution. Allow the disinfectant solution to dry in the air before using the device.*

***NOTE:** water-based disinfecting solutions with a neutral pH are highly recommended; some alcohol-based disinfection solutions may be harmful and discolour or otherwise damage plastic materials.*

8.3 - STERILIZATION PROCEDURE

Warning: carry out sterilization using only a steam autoclave. Do not use other type of sterilization procedure (such as dry, radiation, ethylene oxide, gas, low temperature plasma etc.)

In order to avoid bacterial or viral infection, always clean, disinfect and sterilise the following components after each treatment:

- 1 handpiece (Fig. 2),
- 2 inserts (Fig. 3 - Ref. J),
- 3 torque wrench (Fig. 3 - Ref. C),
- 4 pump tube (Fig. 3 - Ref. P),
- 5 cord/peristaltic pump tube connection (Fig. 3 - Ref. E),
- 6 Handpiece support (Fig. 3 - Ref. O),

The above components are made of materials that will withstand a maximum temperature of 135°C for a maximum of 20 minutes.

8.4 - AUTOCLAVE STERILIZATION OF THE HANDPIECE

Warning: handpiece and cord cannot be detached;

Warning: don't dip the handpiece into disinfectant solutions or liquids of any other kind since this could damage it;

Warning: don't sterilize the handpiece with the insert screwed into it;

Warning: the electric contacts of the connectors of the handpiece and of the cord must be dry;

At the end of the sterilization cycle and before connecting the handpiece to the cord, make sure that the electric contacts of both connectors are completely dry, if necessary, dry the contacts by blowing air on to them with the syringe.

Warning: after completing the sterilization cycle, allow the handpiece to dry completely before using it.

Note: water-based disinfecting solutions with a neutral pH 7 are highly recommended. Some alcohol-based disinfecting solutions may be harmful and discolour to otherwise damage plastic materials.

- 1 Clean the handpiece carefully paying special attention to the threaded pin onto which the inserts are screwed and to the adjacent ring-shaped cavity;
- 2 Disinfect the handpiece using a cloth moistened with a mild disinfectant solution having a neutral pH 7;
- 3 Dry the electric contacts by blowing air onto them with the syringe;
- 4 Seal the handpiece in an individual disposable bag (without any inserts);
- 5 Sterilize the handpiece in the autoclave. Before connecting the handpiece to the main unit, make sure that the electric contacts of both connectors are complete dry, if necessary; dry the contacts by blowing air on to them with the syringe.

8.5 - AUTOCLAVE STERILIZATION OF THE INSERTS

- 1 Clean the insert (preferably in an ultrasonic tank) and rinse it in distilled water;
- 2 Dry the insert;
- 3 Disinfect the insert with a mild disinfectant solution having a neutral pH 7 and dry it carefully.

Warning: before starting the sterilization cycle make sure that the insert is completely dry also inside, to do this, blow air through the internal hole with the syringe.

- 4 Seal the inserts inside individual disposable bags;
- 5 Sterilize in autoclave the inserts.

8.6 - AUTOCLAVE STERILIZATION OF THE TORQUE WRENCH

- 1 Clean the wrench;
- 2 disinfect the wrench with a mild disinfectant solution having a neutral pH 7 and dry it thoroughly;
- 3 seal the wrench inside an individual disposable bag;
- 4 Sterilize in autoclave the wrench.

8.7 - AUTOCLAVE STERILIZATION OF THE PERISTALTIC PUMP TUBE

- 1 Clean the peristaltic pump tube;
- 2 Disinfect with a mild disinfectant solution having a neutral pH 7 and dry it thoroughly;
- 3 Seal the tube inside an individual disposable bag;
- 4 Sterilize in autoclave the tube.

8.8 - AUTOCLAVE STERILIZATION OF THE CONNECTOR BETWEEN CORD AND PERISTALTIC PUMP TUBE CONNECTION

- 1 Clean the connector between cord and peristaltic pump tube connection;
- 2 Disinfect the connector with a mild disinfectant solution having a neutral pH 7 and dry it thoroughly;
- 3 Seal the connector inside an individual disposable bag;
- 4 Sterilize in autoclave the connector.

8.9 - AUTOCLAVE STERILIZATION OF THE HANDPIECE SUPPORT

- 1 Clean the handpiece support;
- 2 Disinfect with a mild disinfectant solution having a neutral pH 7 and dry it thoroughly;
- 3 Seal the connection inside an individual disposable bag;
- 4 Sterilize in autoclave the handpiece support.

9 - REGULAR MAINTENANCE

- 1 Handle this device gently, keep away from the shocks, and install and store away from sunlight.
- 2 Do not store together with poisonous, caustic, explosive or combustible items.
- 3 Store in an environment with humidity under 80%, pressure between 50kPa~106kPa, temperature between -10°C ~ +50°C.
- 4 If the device is not used for a long time, it is better to connect the electricity and water supply one time per month, 5 minutes per time.
- 5 Disconnect the device from the power main when not in use.

Danger: check regularly that the power cable is intact, if it is damaged, replace it with a Carlo De Giorgi spare.

10 - REPLACEMENT OF THE FUSES

Danger: switch off the apparatus.


Fig. 11

Always turn off the device by means of the switch (Fig. 5 - Ref. B) and disconnect it from the power outlet before carrying out the following maintenance activities.

- 1 Insert the flat tip of a screwdriver into the recess in the fuse compartment below the power socket and use it as a lever (Fig.11 - Ref. A).
- 2 Pull out the fuse compartment (Fig.11 - Ref. B).
- 3 **Danger:** Replace the fuses, using fuses of the type indicated in Chapter 15 – Par. 5.
- 4 Put the compartment back into place (Fig.11 - Ref. B).

11 - DISPOSAL PROCEDURES AND PRECAUTIONS

Danger: Hospital waste

Treat the following items as hospital waste:

- Inserts, when worn or broken.
- Irrigator, after each treatment.
- Tube of the peristaltic pump, after 8 sterilizing cycles.
- Torque wrench for tightening inserts, when worn or broken.

12 - INSERTS

12.1 - SHARP INSERTS

The sharp edges of these inserts can be used to treat bone structures efficiently and effectively. Sharp inserts are used in osteotomy and osteoplasty when a fine and well-defined cut in the bone structure concerned is required, there are also inserts with sharp edges for osteoplasty techniques and for removing bone fragments.

12.2 - SMOOTHING INSERTS

The smoothing inserts have surfaces shaped in such a way that they can be used to work the bone structures with precision and in a controlled manner. Smoothing inserts are used in osteotomy when it is necessary to prepare difficult and delicate structures such as those for preparing a maxillary sinus window or to complete preparation of the site of an implant.

12.3 - BLUNT INSERTS

Blunt inserts are used for separating the soft tissues, for example for detaching Schneider's membrane or for lateralizing nerves. In periodontology, these inserts are used to smooth the root surfaces.

13 - SYMBOLS

Trademark.

Caution, consult accompanying documents.

Type B applied part.

Use indoor only.

Alternating current.


Can be autoclaved.


Socket for the foot switch.


CE marked product.

Temperature limitation.

Atmospheric pressure for storage.

Humidity limitation.


Appliance compliance WEEE directive.

Date of manufacture.

Manufacturer.

Authorised Representative in the EUROPEAN COMMUNITY.

ISO 9001:2015

Certificated by international quality control system of ISO.

ISO 13485:2012

Certificated by medical instrument manufacturer quality control system of ISO.

14 - TROUBLESHOOTING

If the device does not seem to be working properly, read the instruction manual again and then check the following table:

Problem	Possible cause	Solution
The device does not turn on when the switch is positioned on ON.	The connector on the end of the power cable is plugged into the socket on the rear of the device incorrectly.	Check that the power cable is firmly connected.
	The power cable is faulty.	Check that the power outlet is working properly. Replace the power cable.
	The fuses blew out.	Replace the fuses.
The device is switched ON but will not work, no errors are shown on the display.	The connector of the footswitch is not properly plugged into the socket.	Insert the footswitch connector properly.
	The footswitch is broken.	Contact the nearest dealer or authorized Carlo De Giorgi service Centre.
A faint whistle can be heard coming from the Ultrasurgery handpiece during operation.	The insert is not correctly tightened onto the scaler handpiece.	Unscrew the insert and screw it back into place correctly.
The device is switched on but does not work, the message WARN appears on the display.	The insert is not fitted correctly into the handpiece.	Unscrew the insert and screw it back into place correctly.
	The insert is worn, broken or deformed.	Replace the insert.
	The connector of the cord is wet.	Dry the connectors.

Problem	Possible cause	Solution
The device is switched on but will not work, the message WARN appears on the display.	Cord not connected to the device.	Connect the cord to the device.
	Lack of continuity of a lead in the cord.	Contact the nearest dealer or authorized Carlo De Giorgi service centre.
	Handpiece failure.	Contact the nearest dealer or authorized Carlo De Giorgi service centre.
	Malfunctioning of the tuning circuit.	Contact the nearest dealer or authorized Carlo De Giorgi service centre.
No liquid comes out of the insert during operation.	The insert is of the type with no through-flow of liquid.	Use an insert of the type with through-flow of liquid.
	The bag of liquid is empty.	Replace the bag with a full one.
	The lid of pump that connected with the water tube is open.	Close the cover.
	The tubes of the drip system and of the pump have not been correctly installed.	Check the connections of the tubes.
	The insert is clogged.	Free the passage in the insert through which the water passes.
	The handpiece is clogged.	Contact the nearest dealer or authorized Carlo De Giorgi service centre.

Problem	Possible cause	Solution
The device is working properly, but the pump is being forced.	Too much pressure by the impeller on the tube in the peristaltic pump.	Check that the tube in the peristaltic pump has been correctly inserted.
The pump is running correctly but when it stops liquid comes out of the handpiece.	The lid of the peristaltic pump is not closed properly.	Make sure that the lid of the peristaltic pump is properly closed.
Insufficient power.	The insert is not correctly fitted to the handpiece (the message WARN appears on the display).	Unscrew the insert and screw it back into place correctly.
	The insert is worn, broken or deformed (the message WARN appears on the display).	Replace the insert.

15 - TECHNICAL DATA

- 1 Device in accordance with Directive 93/42EEC.
- 2 According to EN60529:
 - IPX1 (device)
 - IPX8 (footswitch)
- 3 Device for intermittent operation: Max T_{on} 60s, T_{off} 10s
- 4 Power-supply voltage: 200-240V~ 50/60Hz 360mA
- 5 Fuses: 2×0.6AT 250V
- 6 Working frequency: 24kHz ~ 29.5kHz
- 7 Flow: 25 ~ 100ml/m
- 8 Protection systems and tripping time of the APC:
 - No handpiece connected: 10ms
 - Cord interrupted: 10ms
 - Inserts broken or not correctly tightened: <500ms
 - Protection by discharge to earth: 10ms.
- 9 Warnings: Front display shows WARN XX (see section 7.4 for warnings description of warning and chapter 14 for troubleshooting).
- 10 Operation environment:
 - a) Environment temperature: +10°C ~ +40°C
 - b) Relative humidity: 30% ~ 75%
 - c) Atmosphere pressure: 700 hPa ~ 1060 hPa
 - d) Temperature in the water inlet of water-cooling equipment is not higher than 25°C.
- 11 Delivery and store environment: store in the environment of humidity under 80%, pressure around: 50kPa ~ 106kPa, temperature between -10°C ~ +50°C.
- 12 Handpiece cable: Suggested life cycle - less than 100 sterilization cycles.
- 13 Pump tube: Suggested life cycle - less than 8 sterilization cycles.
- 14 Size of Main unit: 330×254×167 mm.
- 15 Weight of main unit: 3.8kg.
- 16 Type of protection against electric shock: Class I equipment.
- 17 Degree of protection against electric shock: Type B applied part.

16 - AFTER SERVICE

We offer one year, free repair to the equipment according to the warranty card. The repair of the equipment should be carried out by our professional technicians. We are not responsible for any damage caused by the unauthorized non-professional personnel.

17 - ENVIRONMENTAL PROTECTION

There is not any harm factor in our product. You can dispose of it based on local laws.

18 - MANUFACTURER'S RIGHT

We reserve the right to change the design of the equipment, the technical specifications, fittings, the instruction manual and the content of the original packing list at any time without notice. If there are some differences between blueprint and real equipment, take the real equipment as the norm.

19 - PLEASE CONTACT, FOR TECHNICAL DATA

Carlo De Giorgi Srl
Via Tonale 1/3
20021 Baranzate, Milan, Italy
Tel. 02 356 15 43 Fax 02 356 18 08
<http://www.degiorgi.it>

20 - DECLARATION OF CONFORMITY

20.1 - PRODUCT CONFORMITY THE FOLLOWING STANDARDS

EN 60601-1:1990+A1:1993+A2:1995+A13:1996
EN 60601-1-2:2001/EN 60601-1-4:2001
EN 61205:1994/EN ISO 22374:2005
EN ISO 1497:2000+A1:2003
EN 980:2003/ISO 9687:1993/EN 1041:1998
EN ISO 17664:2004/EN ISO 17665-1:2006
EN ISO 10993-1:2003/EN ISO 10993-5:1999/EN ISO
10993-10:2003+A1:2003

20.2 - DECLARATION OF CONFORMITY - EMC

Electromagnetic emissions - Guidance and manufacturer's declaration		
The device Ultrasurgery is intended for use in the electromagnetic environment specified below. The customer or the user of the device Ultrasurgery should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The model Ultrasurgery uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The model Ultrasurgery is suitable for used in domestic establishment and in establishment directly connected to allow voltage power supply network which supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000 - 3 - 0	Not applicable	

Electromagnetic immunity – Declaration & Guidance			
The device Ultrasurgery is intended for use in the electromagnetic environment specified below. The customer or the user of the device Ultrasurgery should assure that is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment- guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 6 kV air	± 6 kV contact ± 6 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for Input/output lines	± 2 kV for power supply lines	Mains power quality should be that of typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 2 kV common mode	Mains power quality should be that of typical commercial or hospital environment.
Voltage dips. Short interruptions and voltage variations on power supply input lines IES 61000-4-11	$< 5\% U_T$ ($>95\%$ dip in U_T) for 0.5 cycle $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 cycles $<5\% U_T$ ($>95\%$ dip in U_T) for 5 sec.	$< 5\% U_T$ ($>95\%$ dip in U_T) for 0.5 cycle $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 cycles $<5\% U_T$ ($>95\%$ dip in U_T) for 5 sec.	Mains power quality should be that of typical commercial or hospital environment. If the user of the model Ultrasurgery requires continued operation during power mains interruptions, it is recommended that the model Ultrasurgery should be connected to an autonomous power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 6100-4-8	3 A/m	Not applicable	Not applicable

NOTE UT is the a. c. mains voltage prior to application of the test level.

Electromagnetic immunity – Declaration & Guidance

The device Ultrasurgery is intended for using in the electromagnetic environment specified below. The customer or the user of the device should make sure that it is used in such an environment.

Immunity Test	Test 60601 IEC level	Compliance level	Guidance - Electromagnetic environment
Conducted R F IEC 61000-4-6	3Vrms 150kHz to 80 MHz	3V	Portable and mobile RF communications equipment should be used no closer to any part of the model Ultrasurgery, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \times p^{1/2}$ $d = 1.2 \times p^{1/2}$ 80 MHz to 800 MHz
Radiated RF IEC 61000-4-3	3V/m 80 MHz to 2.5 GHz	3V/m	$D = 2.3 \times p^{1/2}$ 80 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location is which the model Ultrasurgery is used exceeds the applicable RF compliance level above, the model Ultrasurgery should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the model Ultrasurgery.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m

Recommended separation distances between portable and mobile RF communications equipment and the model Ultrasurgery

The device Ultrasurgery is intend for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the device Ultrasurgery can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device Ultrasurgery as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150kHz to 80 MHz $d = 1.2 \times p^{1/2}$	80 MHz to 800 MHz $d = 1.2 \times p^{1/2}$	80 MHz to 2.5 GHz $D = 2.3 \times p^{1/2}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitter rated at a maximum output power not listed above, the recommended separation distanced in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) accordable to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not be applied in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

The device has been tested and homologated in accordance with EN 60601-1-2 for EMC. This does not guarantee in any way that this device will not be effected by electromagnetic interference.

Avoid using the device in high electromagnetic environment.

21 - GUARANTEE

- 1 Before being placed on the market, all Carlo De Giorgi equipment undergoes a thorough final check to ensure that it is in proper working order.
- 2 Carlo De Giorgi guarantees its products, purchased new from a Carlo De Giorgi dealer or importer, to be free from manufacturing or material defects for:
 - 1 YEAR from the date of purchase for the device;
 - 1 YEAR from the date of purchase for the handpiece with its cord.
- 3 Throughout the warranty period, Carlo De Giorgi undertakes to repair (or, at their sole discretion, to replace) free of charge any parts that, in their opinion, are faulty.
Complete replacement of Carlo De Giorgi products is excluded.
- 4 Carlo De Giorgi cannot accept any liability for direct or incidental damage or personal injury in the following cases:
 - (1) If the equipment is used for purposes other than that for which it is intended.
 - (2) If the equipment is not used in accordance with all the instruction and requirement described in this manual.
 - (3) If the wiring system in the room where the equipment is used does not comply with the applicable standards and appropriate requirements.
 - (4) If any assemble operations, extensions, settings, alterations or repairs have been carried out by personnel not authorised by Carlo De Giorgi.
 - (5) If the environmental conditions in which the device is kept and stored do not comply with the requirements indicated in the chapter on technical specifications.
- 5 Accidental damages due to transport, incorrect use or carelessness or to connection to power supplies other than as envisaged and damage to the signalling lamps handpiece and all accessories are excluded from the warranty.
The warranty will no longer apply if the apparatus has been tampered with or repaired by unauthorised personnel.

6 Warning:

The warranty is valid only if the warranty slip enclosed with the product has been completed in full and returned to us or, if appropriate, to your Carlo De Giorgi dealer or importer within 20 days from the date of purchase, as proven by the consignment note/invoice issued by the dealer/importer.

In order to benefit from the warranty service, the customer must return the apparatus to be repaired to the Carlo De Giorgi dealer/importer from which it was purchased, at his own expense.

7 The apparatus should be returned suitably packed (possibly in its original packing material), accompanied by all the accessories and by the following information:

- a) Owner's details, including his telephone number;
- b) Details of the dealer/importer;
- c) Photocopy of the consignment note/purchase invoice of the apparatus issued to the owner and indicating, in addition to the date, also the name of the apparatus and its serial number;
- d) A description of the problem.

9 Transport and any damages caused during transport are not covered by the warranty.

In the event of failure due to accidents or improper use, or if the warranty has lapsed, repairs to Carlo De Giorgi products will be charged on the basis of the actual cost of the materials and labour required for such repairs.

carlo de giorgi s.r.l.

Via Tonale, 1 – 20021 Baranzate (MI) – Italy
Tel. +39.02.356.15.43 r.a. – Fax +39.02.356.18.08
<http://www.degiorgi.it> E-Mail: info@degiorgi.it