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ENGLISH MANUAL

Thank you for purchasing our product; please keep the instruction leaflet for future reference.

Products' appearance and specifications are subject to change without notice, due to continuous improvement.

WARRANTY

Carlo De Giorgi S.r.l. guarantees the Handpiece, Control box and Foot pedal for 1 year (12 months).

Warranty does not include normal wear of the bearings, collet chuck spindle and armature.

Carlo De Giorgi S.r.l. is not responsible for accidents or malfunctions due to incorrect use of the product.

The users are responsible for the operations and maintenance of the product; in case of malfunctions, please report the model name, serial number, malfunction encountered and contact details for a quicker and more precise service.

VOIDING OF THE WARRANTY

The following occurrences will void the warranty for the product:


- Use with the wrong Voltage input;
- Misuse, against the indications of the instruction manual;
- Damage resulting from falls;
- Use of consumable parts or accessories not approved by Carlo De Giorgi S.r.l.;
- Repair by technicians not authorized by Carlo De Giorgi S.r.l.


PRODUCT SPECIFICATION

PRIME 115/PRIME 407 - CONTROL BOX

INPUT POWER SOURCE	FUSE	WEIGHT	MAXIMUM LOAD CURRENT	DIMENSIONS (mm)
AC 110~120V 50/60 Hz	4A	2.32kg	2A	W. 160 L. 210 H. 100
AC 220~240V 50/60 Hz	2A	2.32kg	1A	


HANDPIECE PRIME 115			
OUTPUT	210W		
SPEED	Max 50,000 RPM		
TORQUE	6.8 N cm		
WEIGHT	195 g		
DIMENSIONS (mm)	L 162	D 20	H 26
HANDPIECE PRIME 407			
OUTPUT	230W		
SPEED	Max 50,000 RPM		
TORQUE	7.8 N cm		
WEIGHT	205 g		
DIMENSIONS (mm)	L 162	D 20	H 28


PRIME 202 - CONTROL BOX

INPUT POWER SOURCE	FUSE	WEIGHT	MAXIMUM LOAD CURRENT	DIMENSIONS (mm)
AC 110~120V 50/60 Hz	4A	2.32kg	1.2A	W. 130 L. 175 H. 77
AC 220~240V 50/60 Hz	2A	2.32kg	0.6A	


HANDPIECE PRIME 202				
OUTPUT	80W			
SPEED	Max 35,000 RPM			
TORQUE	4.8 N cm			
WEIGHT	215 g			
DIMENSIONS (mm)	L 157	D 16	H 27	


COMPONENTS LIST

PRIME 115/407

- 1) Control box;
- 2) RPM display;
- 3) Volume dial;
- 4) Motor ON/OFF button;
- 5) Hand/foot button;
- 6) Forward/reverse button;
- 7) Motor connector;
- 8) Power switch;
- 9) Foot pedal connector;
- 10) Fuse box;
- 11) Power connection;
- 12) Power cord;
- 13) Foot speed control pedal;
- 14) Handpiece stand;
- 15) Handpiece;
- 16) Wrench.


PRIME 202

- 1) Control box;
- 2) Volume dial;
- 3) Motor connector;
- 4) Forward/reverse button;
- 5) Power switch;
- 6) Fuse box;
- 7) Power connection;
- 8) Foot pedal connector;
- 9) Power cord;
- 10) Foot speed control pedal;
- 11) Handpiece stand;
- 12) Wrench;
- 13) Handpiece;
- 14) Power light.

SETUP AND OPERATION – PRIME 115/407

WARNING

To avoid malfunctions, please make sure that the power switch is off before connecting the motor and foot pedal.

- 1) Verify that the power switch is off.
Connect the power cord (A).
Connect the foot control pedal to the foot pedal connector (B).
Connect the handpiece to the handpiece connector (C).
- 2) Put the speed control dial on the control box at "0".
- 3) Check the proper lock of the collet chuck, by turning the handpiece handle or pulling out the bur.
- 4) Hold the handpiece with one hand and turn on the power switch of the control box with the other hand.
The power switch lights green, and the RPM display turns on.
- 5) To power off the device, turn off the power switch.
The green light on the power switch, the button lights and the RPM display turn off.

- 6) Cycle between hand and foot control by pressing the “Hand/Foot” select button.

- 7) Cycle between left and right rotation by pressing the “Right/Left” select button.

- 8) Press the Motor ON/OFF button to start or stop the motor.

- 9) Hand control icon.

Foot control Icon.

Verify the correct functioning of the handpiece while adjusting the speed dial of the control box after pushing the Motor ON/OFF button.


The Motor speed can be changed up to the limit preset by the speed dial of the control box. The speed will increase with the pressure on the pedal.

IN CASE OF OVERLOAD

In case of overload, after 5 seconds, the control box stops and shows the error message (---) with beeper sound. To restart, push the Motor button if in hand control mode, or take off the foot from the pedal if in foot control mode.

SETUP AND OPERATION – PRIME 202

WARNING

To avoid malfunctions, please make sure that the power switch is off before connecting the motor and foot pedal.

- 1) Verify that the power switch is off.

Connect the power cord (A).

Connect the foot control pedal to the foot pedal connector (B).

Connect the handpiece to the handpiece connector (C).

- 2) Put the speed control dial on the control box at "0".
- 3) Check the proper lock of the collet chuck, by turning the handpiece handle or pulling out the bur.
- 4) Hold the handpiece with one hand and turn on the power switch of the control box with the other hand.
The red light on the display turns on.
- 5) To power off the device, turn off the power switch.
The red light on the display turns off.

- 6) Change between left and right rotation by pressing the “Right/Left” switch, on the side of the control box.
- 7) When the pedal is inserted, the foot control is active. The Motor speed can be changed up to the limit preset by the speed dial of the control box. The speed will increase with the pressure on the pedal.

When the pedal is removed, the hand control is active. The speed will increase with the rotation of the speed dial.

IN CASE OF OVERLOAD

In case of overload, after 5 seconds, the control box stops. To restart, push the Motor button if in hand control mode, or take off the foot from the pedal if in foot control mode.

MAINTENANCE ADVICE

- 1) Keep the handpiece body clean.
- 2) Remove any dust on the handpiece with smooth cloth soaked with isopropyl alcohol, being careful not to let any liquid inside the handpiece.
 - Do not apply oil or water to the device for cleaning, as it may damage the bearings or rust the inner parts.
- 3) Air cleaning is possible, but strong air pressure may damage the handpiece.


ISSUE	CHECK	SOLUTION
RPM display doesn't appear	Check if the control box is well connected with the power outlet.	Check if the power cord connects well from the outlet to connector.
	Check if the power is on and the light of the switch is on.	Change the fuse inside (refer to the fuse change section in the repair manual). Change the fuse wire in the inlet if it is cut.
Buttons are not working	Check if the device does not work after the beeping sound.	Send to the distributor or manufacturer for repairs.
Handpiece heats up	The handpiece is warm.	The handpiece heats up slightly when it is used continuously. Suspend the operations for a short break.
	The handpiece is very hot.	The handpiece heats up considerably if used close to overload conditions. If the handpiece heats up considerably, and the overload condition does not occur, send to the distributor or manufacturer for repairs.
Handpiece doesn't rotate	Check if the chuck is open.	Close the chuck.
	Check if handpiece is connected correctly to the control box.	If the handpiece is connected correctly but it does not work, send to the distributor or manufacturer for repairs.

CAUTION

- 1) The maximum working RPM for the device varies according to the bur diameter.

Please keep the RPM under the limits described in the following diagram to avoid malfunctions.


- 2) If any problem occurs in the handpiece, do not use it until repaired.
- 3) Make sure that the switch is off before putting the power plug in the socket, as the device might start unexpectedly otherwise.
- 4) A circuit breaker is provided to protect the motor. The breaker activates when the load on the motor is too high. Please use the product within allowable ranges.

- 5) Do not open the chuck while operating.
- 6) Keep a bur in the handpiece collet chuck while not in use.
- 7) Do not put the handpiece down before it stops completely.
- 8) Be sure to remove the power plug from the power outlet before detaching the control box.
- 9) Use the specified voltage.
- 10) Use accessories suitable for this device.
- 11) Wear proper protective clothes and equipment, including a dust mask and safety goggles. Use a proper suction unit and a ventilator.
- 12) Keep the workplace clean. Debris may cause an accident.
- 13) The device keeps the last used direction and hand/foot selections. Keep this in mind when turning the device on.
- 14) Do not use the control box in an excessively humid environment.
- 15) Check for any issues, including vibrations of the bur, at low speeds, before turning on high speeds.
- 16) Do not turn on the power of the control box while the collet chuck of the handpiece is open.


PRIME 407/115 – BEARING CHANGE








PRIME 202 – BEARING CHANGE

BEARINGS 1, 2


BEARINGS 3, 4


COLLET CHUCK CHANGE

PRIME 115/407

PRIME 202

CE DECLARATION OF CONFORMITY

The undersigned:

CARLO DE GIORGI S.r.l.

Via Tonale, 1

20021 BARANZATE (MI) - I

declares, under its own responsibility, that the apparatus:

Type: **LABORATORY MICROMOTOR PRIME 115/407/202**

described as:

laboratory micromotor for roughing, polishing and finishing materials utilized in orthodontic laboratories, jewellers workshops and tool shops,

conforms to the Laws implementing the requirements of the "Machinery Directive" 2006/42/CE, the "Low Voltage Directive" 2014/35/UE and the "Electromagnetic Compatibility" Directive 2014/30/UE, in accordance with the applicable technical standards EN ISO 12100, EN 61010-1, EN 61029-1, EN 61029-2-4, EN 60326-1, EN 55014-1, EN 55014-2.

Signed:

De Giorgi Ariberto

CEO

carlo de giorgi

Via Tonale 1

20021 Baranzate - Milano

Tel 02 3561543

Fax 02 3561808

Web site: www.degiorgi.it

Mail: info@degiorgi.it