“ZEUS ”

FORNO PER PRE-RISCALDO

INSTALLATION USE MAINTENANCE

CARLO DE GIORGI S.r.l.

Via Tonale, 1

20021 Baranzate - MI - Italy

Phone: +39 02 3561543 - Fax: +39 02 3561808

E-mail: info@degiorgi.it

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1.1 - Scopo del manuale

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1 - GENERAL INFORMATION

1.1 Aim of the manual

This manual has been written by the manufacturer and it is an integral part of the machine.


Its information are addressed to the user and contain all Safety Recommendations.

Before using this machine, it is recommended that you read this manual carefully, especially if it is the first time, to ensure you know the control devices, their function and position. It is also advisable to make performance tests.

This manual must be kept for future references.


1.2 Identification of the machine manufacturer

On the machine you can find the identification label of the manufacturer and some information regarding technical data as shown in Pic. 1.

1.3 – Technical information

The cylinder pre-heating oven “ZEUS” is managed by a microprocessor electronic control and it has been designed for the firing and pre-heating of the coating of the cylinders to cast metals in dental technicians and goldsmiths labs. The oven can memorise up to 7 heating cycles (or programs) which are kept in the memory until the operator decides to modify them. Each program contains 4 phases and for each phase it is possible to set 3 variables: temperature, increase rate, holding time. It is also possible to set an Ignition delay (up to 100 hours), fumes aspiration fan Off temperature (fume chimney) and a final temperature holding time at the end of the cycle.

Legenda (PIC. 2):

A	Muffle	E	Feeder 230V 50Hz
B	Electronic Control	F	Fumes Aspiration Chimney (Optional)
C	Door		
D	Refractory Chamber		

1.4 - Safety devices

The parts of the machine supplied by electricity are protected by fixed guards to avoid access by the operator. Only specialized and authorized personell must have access to these parts for extraordinary maintainance and repairs.

1.5 - Technical Data

	ZEUS S	ZEUS M	ZEUS L
Power Supply	230V-50/60Hz	230V - 50/60Hz	230V - 50/60Hz
Resistances Power	1600W	2000W	2500W
Chamber dimensions	cm15x15x10	cm18x23x11.5	cm23x30x15
External DimensionsLxPxH	cm32x40x50	cm36x46x52	cm42x53x55
Weight – Kg	25	35	45
Number of Programmes:	7	7	7
Number of Phases	4	4	4
Ignition Delay	Si	Si	Si
Fuses	16A	16A	16A

2 - MOVING AND INSTALLING

2.1 - Packaging and unpackaging

The packaging is as follows:

- Polystyrene plates K10 protecting the machine on sides, top and bottom.
- Hard cardboard cover
- Strap

2.2 Loading and unloading

This machine must be moved by carts or manually by at least 2 people.

While moving the machine avoid absolutely any kind of bumping, dropping or tilting: they could seriously damage it. Obviously, the manufacturer is not responsible for damages caused by droppings, improper use and maintenance which are not in strict accordance with the manufacturer's instructions contained in this manual.

Packaging waste disposal must be done with respect to the environment and the law in force.

2.3 Installation instructions

The machine must be located on a safe place and on a horizontal position. Make sure there is sufficient air where you position the machine.

It is duty of the user to ensure that the electric plant is in accordance with the safety laws in force. It is particular important to make sure that the grounding is well efficient. Furthermore, it is important to verify the network voltage: in case the voltage is too low (lower than 210V), the engine will start too slow and it might be necessary to install a voltage stabilizer.

After having placed the oven and gone through the above due verifications, follow these instrutions:


- The general switch must be in the OFF (0) position.
- Connect the power supply plug to a 220-230 V AC outlet.

3 - USE INFORMATION

3.1 – Performance and Use

The oven as said earlier can memorise up to 7 pre-heating cycles (or programs) which are kept in the memory until the operator decides to modify them. Each program contains 4 phases and for each phase it is possible to set 3 variables: temperature, increase rate, holding time. It is also possible to set an Ignition delay and the fumes aspiration fan Off temperature (fume chimney). Furthermore it is possible to set a final temperature holding time at the end of the cycle: this is to keep the cylinders at a certain temperature and allow some time for the operator to be ready with the next work.

3.1 - Panel - Operating description (Pic.3)

PIC.3

1- RISE TEMPERATURE LED: is lit during the heating phases.

2- RISE RAMP LED: is lit when you set the programming of the rise ramp and during the heating.

3- STAY AND HOLDING LED: is lit during the ignition delay and during the temperature holding phases.

4- DISPLAY: when programming it indicates the values of the different functions that are modified (ignition delay, phase temperature, rise speed, holding time). During the working process it indicates the values of the different functions of the programmed phases (i.e. during the heating it indicates the temperature it has to arrive to, during the holding it indicates the set time with a value that decreases as time elapses).

5- DISPLAY: when programming it indicates the phase you are in - during working process it indicates the inner temperature of the chamber.

6- In this area are indicated the different phases of each programme. Both during programming and working processes, as you change from one phase to the next the related led will light up.

7- In this area the seven memorized programmes are visualised. The programme you are working on is indicated by its lit led.

8- LED WORK: it is lit when the oven is working.

9- LED STOP: it is lit when the oven is not working.

10- START/STOP BUTTON: by pushing it you turn on the oven (which means that an heating cycle will begin); by pushing it when the oven is on, you turn it off.

11- ENTER: when programming it confirms the set values.

12- SET/EXIT BUTTON: by pushing it twice you enter the programming phase, by pushing it again you exit that phase. If you have modified the visualised values in order to exit you have to push first ENTER and then EXIT.

13- SELECTION BUTTONS: you use it to select the programmes, the phases and to modify the set values.

14- PHASE BUTTON: you can choose a different phase from the one working. Use it like the button 15 Prog (push it and then select it with 13).

15- PROG BUTTON: you can go from one programme to another. In fact, if you push it the programme working in that moment will start flashing, then by pushing arrow 13 you can change to a different programme.

3.3 -Working Cycle - Functions

The oven "ZEUS" is managed by a microprocessor electronic control and it has been designed to perform heating or cooling cycles or "programmes" chosen and set by the technician. For each programme it is possible to set the following:

- Ignition Delay (in hours and minutes)
- N° 4 phases (Max) of heating, for each phase it is possible to set 3 variables: Temperature, Increase Rate (in C°/min) and Holding Time Temperature. It is also possible to set cycles with a number of phases lower than 4. In order to eliminate a phase it is sufficient to eliminate the 3 values regarding the phase itself;
- Fumes aspiration fan Off temperature;
- Final temperature holding time at the end of the cycle.

Programming


Choose the programme to modify or set by pushing the button 15 PROG: one of the 7 programmes leds will light up.

With the button 13 move from 1 to 7 or viceversa until you reach the desired programme.


Push the button 12 SET/EXIT twice.

After having pushed it once the lettering PAR will be displayed.

After having pushed it twice the following lettering will be displayed:


You can set an ignition delay by using the arrows 13 (in hours and minutes), then confirm with ENTER passing to phase 1:
The following will be displayed:

Set the temperature of Phase 1 with the arrows 13 and confirm with ENTER. The following will be displayed:


Set the increase rate speed (in C°/min) with the arrows 13 and confirm with ENTER. The following will be displayed:


Set the holding time (in hours/minutes) with the arrows 13 and confirm with ENTER; you enter now phase 2. Continue like this until phase 4.

Once phase 4 has been set (that is after having pushed ENTER for the holding time of phase 4) the following will be displayed:

At this point you are asked to set the fumes aspiration fan Off temperature (this is valid only if the oven is supplied with a chemney and the fan starts at the beginning of the heating): set now, with the arrows 13, the value of the temperature you want the fan to turn off, then confirm with ENTER and the following will be displayed:

Now you are asked to set the holding time: that is the time that the temperature inside the chamber must be maintained at the end of the programme.

During this time the oven will produce an acoustic sound at regular intervals to warn the user that the oven is ready for the casting.

Such time is expressed in hours and minutes, reaching a maximum of 99 hours and 59 minutes.

During the working phase, once the selected time has elapsed, the oven will turn itself off.

START/STOP

To start the oven please follow these steps:

- Choose the programme: push the button 15 PROG and with the arrows 13 find the desired programme.
- If you want an ignition delay, you can set it as follows: calculate the estimated length of the programme, then calculate the time to be set; push the button SET twice, set the delay time, push ENTER and then push again button 13 SET/EXIT to exit the programme.
- Make sure that the door of the oven is safely closed;
- Push the button 10 START/STOP to start the heating cycle;
- Once the firing is completed and the cylinders are removed from the oven, push the button 10 START/STOP to turn off the oven (it is advisable not to let the oven turn off by itself at the end of the holding time since it would be just a waste of energy).

IMPORTANT OBSERVATIONS

If you wish to programme a cycle with a number of phases lower than 4, for example 2, it is always advisable to use the first two phases, that is phase 1 and phase 1 and make ineffective phases 3 and 4: meaning that, for phases 3 and 4, you set the phase temperature equal to the one of phase 2, the increase rate speed equal to phase 2 and the holding time at 0.

By doing this the phase 3 and 4 will be ignored.


During the working phase the display 5 always indicates the real temperature inside the chamber.

Instead, the display 4 indicates during the heating the phase temperature that must be reached, and during holding time, the holding decreasing time left.


SPARE PARTS LIST OVEN ZEUS

Pic. 5- 6

POS	NAME	CODE	ZEUS 1	ZEUS 2	ZEUS 3
1	Cable lead	018.001			
2	Feeder cable	018.002			
3	Fuse holder	018.003			
4	General switch	018.004			
5	Electronic power card	018.005			
6	Micro - Switch	018.006			
7	Panel & electronic card	018.007			
8	Spring Micro	018.008			
9	Handle	018.009			
10	Chimney tube	018.010			
11	Thermocouple	018.011			
12	Flap door spring	018.012			
13	Flap door refractory		018.013S	018.013M	018.013L
14	Frontal refractory		018.014S	018.014M	018.014L
15	Muffle refractory (without resistance)		018.015S	018.015M	018.015L
15	Complete muffle (with resistance)		018.015RS	018.015RM	018.015RL
	Resistance		018.RS	018.RM	018.RL
16	External Chemney (without exhauster)	018.016			
17	Exhauster	018.017			
18	Dewaxing plate		018.018GZS	018.018GZM	018.018GZL


PIC. 5

PIC.6

ELECTRIC PLANT

SHEET TO ORDER SPARE PARTS - OVEN ZEUS

With reference to Pictures N° 5-6 of Manual of Use and Maintenance

Pos	Denomination	Code	Number
1	Cable lead	018.001	
2	Feeder cable	018.002	
3	Fuse holder	018.003	
4	General switch	018.004	
5	Electronic power card	018.005	
6	Micro-Switch	018.006	
7	Panel & electronic card	018.007	
8	Spring Micro	018.008	
9	Handle	018.009	
10	Chimney tube	018.010	
11	Thermocouple	018.011	
12	Flap door spring	018.012	
13	Flap door refractory	018.013 S	
13	Flap door refractory	018.013 M	
13	Flap door refractory	018.013 L	
14	Frontal refractory	018.014 S	
14	Frontal refractory	018.014 M	
14	Frontal refractory	018.014 L	
15	Muffle refractory (without resistance)	018.015 S	
15	Muffle refractory (without resistance)	018.015 M	
15	Muffle refractory (without resistance)	018.015 L	
15	Complete muffle (with resistance)	018.015 RS	
15	Complete muffle (with resistance)	018.015 RM	
15	Complete muffle (with resistance)	018.015 RL	
15	Resistance	018.RS	
15	Resistance	018.RM	
15	Resistance	018.RL	
16	External Chemney	018.16	
17	Exhauster	018.17	
18	Dewaxing plate	018.018GZS	
18	Dewaxing plate	018.018GZM	
18	Dewaxing plate	018.018GZL	