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
declares, under its own responsibility, that the apparatus:

Type: **LABORATORY MICROMOTOR**

described as:

laboratory micromotor for roughing, polishing and finishing materials utilized in orthodontic laboratories, jewellers workshops, tool shops, aesthetic and podiatric centers,

conforms to the Laws implementing the requirements of the "Machinery Directive" 2006/42/CE, the "Low Voltage Directive" 2014/35/UE and the "Electromagnetic Compatibility" Directive 2014/30/UE, in accordance with the applicable technical standards EN ISO 12100, EN 61010-1, EN 61029-1, EN 61029-2-4, EN 60326-1, EN 55014-1, EN 55014-2.

Signed:

De Giorgi Alberto
CEO

Date: 10/10/2016

STANDARDS AND GENERAL PRECAUTIONS

INTRODUCTION

This manual is the property of CARLO DE GIORGI S.r.l. The contents of this manual cannot be reproduced or disclosed to other parties. All rights reserved.

DESCRIPTION OF THE DEVICE

The device is a laboratory micromotor designed for roughing, polishing and finishing of the materials utilized in orthodontic laboratories, jewellers workshops and tool shops.

CAUTIONS BEFORE INSTALLATION

Please note the following important points:

- Do not install in a location where water might be splashed on the product or in a location exposed to such excessive moisture that condensation forms on the product surface.
- Do not install in a location where fire, gas or oxygen are released.

The device must be installed in a dry, well-lit environment in accordance with the requirements of established legislation.

NOTE: The expression "established legislation" refers to the laws in force in the user's country.

REQUESTS FOR SERVICE AND SPARE PARTS

For all electrical maintenance requirements, contact CARLO DE GIORGI at Baranzate, Via Tonale, 1, Milan, Italy.

In case of any faults, notify CARLO DE GIORGI's technical personnel of the exact nature of the problem, specifying also the lot number shown on the label affixed to the bottom of the appliance.

TECHNICAL FEATURES

TECHNICAL FEATURES AND PERFORMANCE

Performance: intermittent duty.

INPUT POWER SOURCE	FUSE	WEIGHT	MAXIMUM LOAD CURRENT	DIMENSIONS (mm)
AC 220~240V 50/60 Hz	1A	0,97kg	1A	W. 101 L. 104 H. 85

	Handpiece 436/00 e 437/00	Handpiece 438/00R
POWER	65W	65W
SPEED	Max 35,000 RPM	Max 35,000 RPM
TORQUE	3.0 N cm	3.0 N cm
WEIGHT	200 g	140 g
DIMENSIONS (mm)	L 159 D 18.9 H 28.2	L 82 D 25

POSITIONING


The appliance is supplied assembled and ready for use after connecting the control devices and the micromotor to the power supply.

The appliance can be moved by hand without requiring lifting devices.

ATTENTION:

Plug into the correct power socket, after having verified the information on the label situated under the device (230V).

COMPONENTS – 436/00**Figure 1 – Micromotor 436/00 components**

- | | |
|---------------------------|-------------------------|
| 1. Control box | 6. Motor connector |
| 2. Handpiece | 7. Speed control dial |
| 3. On/off foot pedal | 8. Foot pedal connector |
| 4. Power switch | 9. Power connector |
| 5. Forward-reverse switch | 10. Fuse holder |

COMPONENTS – 437/00



Figure 2 – Micromotor 437/00 components

- | | |
|-----------------------------|-------------------------|
| 1. Control box | 7. Speed control dial |
| 2. Handpiece | 8. Foot pedal connector |
| 3. Speed control foot pedal | 9. Power connector |
| 4. Power switch | 10. Fuse holder |
| 5. Forward-reverse switch | 11. Pedal/hand selector |
| 6. Motor connector | |

COMPONENTS – 438/00R**Figure 3 – Micromotor 438/00R components**

- | | |
|-------------------------------|-------------------------|
| 1. Control box | 6. Motor connector |
| 2. Motor with intra connector | 7. Speed control dial |
| 3. On/off foot pedal | 8. Foot pedal connector |
| 4. Power switch | 9. Power connector |
| 5. Forward-reverse switch | 10. Fuse holder |

CONNECTION AND PLACING INTO SERVICE

With reference to Figure 1, after positioning the device and checking that the main power switch (4) is set to OFF, connect the following parts:

- foot pedal (3) to connector (8) on the rear of the appliance;
- micromotor (2) to connector (6) on the front of the appliance;
- mains power cable, designed for connection to a mains socket outlet with overcurrent protection devices.

CAUTION:

Before connecting the foot-pedal (3) disconnect the device from the source of power, disconnecting the power cable from the power socket, and connect it again after the foot-pedal (3) is connected to the device (8).

- N.B.** DURING THE INITIAL PERIOD OF USE THE HANDPIECE MAY TEND TO BECOME HOT DURING OPERATION. THIS PHENOMENON DEPENDS ON THE INITIAL ADJUSTMENT OF THE TRANSMISSION ELEMENTS AND WILL DISAPPEAR WITHIN ABOUT ONE WEEK OF NORMAL USE. IF THE PROBLEM SHOULD PERSIST, CONSULT CARLO DE GIORGI.

OPERATION

Keep the device on a table in front of the user.

- 1 - Insert the power cord in the power socket after turning off the power unit switch (4).
- 2 - Connect the handpiece (2) to the output socket (6) of the power unit.
- 3 - Set the RPM to "0" or a low position.
- 4 - Verify that the handpiece chuck is closed by turning the chuck level or pulling out the bar. With the 438/00R model, make sure the chosen handpiece is well connected to the motor, and the chosen tool is well serrated in the handpiece.
- 5 - Turn on the main unit by putting the power unit switch (4) in the "ON" position.
- 6 - Set the RPM with the potentiometer knob (7), according to the machining to be performed.
- 7 - To reverse the direction, turn off the main unit turning off the power unit switch (4) and switch the direction (5).
- 8 - To use the foot pedal (3) connect it to the socket (8) in the rear of the power unit. With the 437/00 model, use the selector (11) on the rear of the power unit.

IMPORTANT:

To use the foot switch (3), set the power unit switch on "OFF" (4).

If the power unit switch (4) is positioned on "ON", the pedal will not work correctly.

INTENDED USE

The laboratory micromotor is designed for use in orthodontic laboratories, jewellers workshops and tool shops, for machining various materials including resins, ceramics, gold, chromium cobalt, etc.


The device is not designed for use in applications that fall within the definitions of the Medical Equipment Directive 93/42/EEC, or any uses other than those specified by the manufacturer.

MAXIMUM WORKING SPEED

The maximum working RPM for the device varies according to the bur diameter.

Please keep the RPM under the limits described in the following diagram to avoid malfunctions.

OPERATOR INSTRUCTIONS

CONTROLS AND POWER SUPPLY UNIT

The laboratory micromotor is equipped with the following command and control devices located on the power supply unit:

- a) main power switch with green indicator (4);
- b) black selector for the tool rotation direction (reverse) (5);
- c) tool rotation footpedal (3);
- d) tool speed slider (7);
- e) power circuit fuse (10).

DIAGNOSTICS

If the appliance is not working first check the condition of the fuse in the rear (10).

If necessary, install a replacement fuse rated as specified to pag. 4.

Before checking the fuse, set the main switch (4) to "OFF" and disconnect the device from the electrical power supply; now check the fuse and replace if necessary.

After the screwcap has been refitted, the power cable can be reconnected.

If malfunctions and faulty operation are persistent, set the main power switch to OFF, unplug the device and consult CARLO DE GIORGI.

When reporting a malfunction, please make sure to explain the exact nature of the problem to our technical personnel.

CAUTION:

Any intervention on the device must be performed on the device isolated from every source of energy.





MAINTENANCE ADVICE





- 1) Keep the handpiece body clean.
- 2) Remove any dust on the handpiece with smooth cloth soaked with isopropyl alcohol, being careful not to let any liquid inside the handpiece.
 - Do not apply oil or water to the device for cleaning, as it may damage the bearings or rust the inner parts.
- 3) Air cleaning is possible, but strong air pressure may damage the handpiece.

ISSUE	CHECK	SOLUTION
RPM display doesn't appear	Check if the control box is well connected with the power outlet.	Check if the power cord connects well from the outlet to connector.
	Check if the power is on and the light of the switch is on.	Change the fuse inside (refer to the fuse change section in the repair manual). Change the fuse wire in the inlet if it is cut.
Buttons are not working	Check if the device does not work after the beeping sound.	Send to the distributor or manufacturer for repairs.
Handpiece heats up	The handpiece is warm.	The handpiece heats up slightly when it is used continuously. Suspend the operations for a short break.
	The handpiece is very hot.	The handpiece heats up considerably if used close to overload conditions. If the handpiece heats up considerably, and the overload condition does not occur, send to the distributor or manufacturer for repairs.
Handpiece doesn't rotate	Check if the chuck is open.	Close the chuck.
	Check if handpiece is connected correctly to the control box.	If the handpiece is connected correctly but it does not work, send to the distributor or manufacturer for repairs.

VISUAL GUIDE FOR BEARING CHANGE – 436/00 AND 437/00


VISUAL GUIDE FOR COLLET CHUCK CHANGE – 436/00 AND 437/00


HANDPIECE EXPLODED DIAGRAM – 436/00 AND 437/00

MOTOR EXPLODED DIAGRAM – 438/00R

DISPOSAL

The crossed out wheeled bin symbol indicates that the apparatus must not be disposed of with other ordinary waste.

The differentiated disposal of this apparatus at the end of its life-cycle is organized and managed by the manufacturer.

When the user wishes to throw this device away he must contact and follow the system adopted by the manufacturer to dispose of it correctly.

Differentiated disposal to enable possible recycling or environmentally compatible elimination helps to limit undesirable effects on health and the environment and promotes the reuse and/or recycling of the materials that make up the appliance.

Legal action may be taken the holder if the apparatus is disposed of illegally, in accordance with the relative legislation in force in the country of use.

carlo de giorgi s.r.l.

Via Tonale, 1 – 20021 Baranzate (MI) – Italy
Tel. +39.02.356.15.43 r.a. – Fax +39.02.356.18.08
<http://www.degiorgi.it> E-Mail: info@degiorgi.it