
PIN 2000

FORAGESSI

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DRILLING MACHINE

USE AND MAINTENANCE

CARLO DE GIORGI S.r.l.

Via Tonale, 1

20021 Baranzate di Bollate - MI - Italy

Phone: +39 02 3561543 - Fax: +39 02 3561808

E-mail: info@degorgi.it

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1.1 - Scopo del manuale

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1 - GENERAL INFORMATION

Aim of the manual

This manual has been written by the manufacturer and it is an integral part of the machine.

Its information are addressed to the user and contain all Safety Recommendations.

Before using this oven, it is recommended that you read this manual carefully, especially the first time, to ensure you know the control devices, their function and position. It is also advisable to make use tests.

This manual must be kept for future references.

1.2 - Identification of the machine manufacturer

On the right of the machine find a CE Mark and the technical.

1.3 - Use

The PIN 2000 Drilling machine is capable of drilling plaster models or Plexiglass and PVC bases or other plastic materials in dental laboratories and in conformity with the safety requirements. The machine is supplied with a Laser pointer system that allows a fast and safe alignment of the perforation.

1.4 - Safety devices

The machine is protected by fixed guards to avoid access to the electrical and mechanical parts by the operator. Only specialized and authorized personell must have access to these parts for extraordinary maintainance and repairs.

1.5 - Technical Data

Dimensions	cmt25L x 26P x 40H
Weight	Kg 15
Power Supply	220 - 230 V
Power	150 W
Rotating Speed	2800 r/min
Tip Diameter	3 mm
Pointer	Laser
Fuses	5A

2 - MOVING AND INSTALLING

2.1 - PACKAGING AND UNPACKAGING

The packaging is as follows:

- Polystyrene plates K10 protecting the machine on sides, top and bottom.
- Hard cardboard cover

2.2 - Loading and unloading

This machine can be moved either by carts or manually.

While moving the machine avoid absolutely any kind of bumping, dropping or tilting: they could seriously damage it. Obviously, the manufacturer is not responsible for damages caused by droppings, improper use and maintenance which are not in strict accordance with the manufacturer's instructions contained in this manual.

Packaging waste disposal must be done with respect to the environment and the law in force.

2.3 - Installation instructions

The machine must be located on a safe place and on a horizontal position. Make sure there is sufficient air where you position the machine.

It is duty of the user to ensure that the electric plant is in accordance with the safety laws in force. It is particular important to make sure that the grounding is well efficient. Furthermore, it is important to verify the network voltage: in case the voltage is too low (lower than 210V), the engine will start too slow and it might be necessary to install a voltage stabilizer. After having placed the machine and gone through the above due verifications, follow these instrutions:

- The general switch must be in the OFF (0) position.
- Connect the power supply plug to a 220-230 V AC outlet, using the provided cable

3 - USE INSTRUCTIONS

3.1 - Operating description (Pic.2)

1 - Cylinder holder with 5A fuses.

2 - Feeder cable 220 V AC

3 - General Switch : by switching it the machine will turn ON and the Display 8 will light up.

4 - Laser ON switch

5 - Button START/STOP : motor ignition

6 - Ignition led: is lit when the general switch is on

7 - Engine functioning led: it flashes when the engine is working

8 - Holding-model tray.

9 - Moving working area: the drilling of the model is done by lowering the working surface with both hands.

10-Regulating Nut: the perforation depth can be regulated by rotating the nut clearly located on the front of the machine.

3.2 - Performance and Use

- Turn ON the machine by using the general switch 3.
- Lay the model or the Plexiglass or the PVC base on the working surface next to the pointer delivery side.
- Push and let the surface 9 go all the way down and check the depth of the milling hole: regulate such depth by rotating the regulating nut 10. By rotating it to the right you'll increase the depth; by rotating it to the left you'll decrease it. Once you have performed this regulation, let the surface free and it will go back to its initial and high position.
- Turn On the motor by using the START/STOP 5 Switch.
- Turn On the Laser pointer by using the Switch 4 .
- Move the model on the surface until the Laser pointer identifies exactly the position of the perforation.
- By holding firmly the model with both hands, push the surface and let it go all the way down: the perforation is obtained in this way.
- Release the pressure on the surface and let it go back to its initial and high position. Proceed then with another perforation.

3.3 - USEFUL RECOMMENDATIONS

AWARE: When the motor is ON, absolutely avoid putting your hands and pushing with them on the working surface: it could cause injuries to your hands!!!

Absolutely avoid regulating the perforation depth when the motor is moving.

When you want to substitute the mill, do as follows: take out the holding-model tray 8, keep in place the spindle with a curved n.10 key, loosen the bush with a n.14 key, take out the old mill and replace it with a new one.

To get rid of chips and dust due to the working process do as follows: take out the holding-model tray 8, suck in with a vacuum cleaner the chips and dust from the inner chamber.

4 - MAINTENANCE INFORMATION

4.1 - Ordinary maintenance

It is advisable to keep the machine free from chips and dust: this easy maintenance operation allows better results and a better machine performance.

Never use any solvents for the cleaning of the machine, it is enough to vacuum the dust.

In order to better clean and suck the dust from the inside of the machine, please do as follows: take the moving working surface 9 apart. To do this (Pic.3) you need to: take the holding-model tray out, loosen of 3 complete turns the screws A, shown in the picture using the supplied key, lift the surface and let it come off its slides, move the surface letting it rotate around the neck of the machine. To put it back in place follow the same instructions backwards.

4.2 - Extraordinary maintenance

To repair the machine or replace some of its parts, the user should consult directly the manufacturer or only specialized technicians.

Absolutely avoid opening any of the fixed guards without proper precautions.

Disconnect the power supply plug from the outlet before any extraordinary maintenance operations.

4.3 - Precautions

Inside the machine there are components subjects to 220V network voltage.

Before performing any operation inside the machine make sure you have turned off the general switch 3 and unplugged the power supply 2 (same precautions must be taken when replacing the fuses).

It is our most concern to point out what has been said in 3.3: the possibility of injuries to the hands exists due to the fact that the engine is always functioning during the working drilling phase. It is fundamental, therefore from the user's end, the maximum attention and the complete observations of the safety regulations.

If the user doesn't strictly follow the rules contained in this manual the manufacturer is lifted from any guarantee and responsibility.

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