*Manuale d'Istruzioni*  
*Instruction's Manual*  
*Manual de Instrucciones*  
*Notice d'Emploi*  
*Bedienungsanleitung*  
*Manual de Instruções*

# RF Ablator

## RadioFrequency

beautyled

### IMPORTANTE

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### IMPORTANT

These operating instructions form an integral part of the equipment and must be available to the operating personnel at all times. All the safety instructions and advice notes are to be observed. Be sure that these operating instructions are furnished together with the equipment when this is transferred to other operating people.

In case of necessity of technical, or other type, assistance contact your own retailer.

*No part of this document could be photocopied, reproduced or translated in other language without the written consent of LED SpA. All rights reserved.*

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| Guidance and manufacturer's declaration – electromagnetic emissions                                                                                                                                   |            |                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The <b>RF ABLATOR</b> is intended for use in the electromagnetic environment specified below.<br>The customer or the user of the <b>RF ABLATOR</b> should assure that is used in such an environment. |            |                                                                                                                                                                                                                                                                                                                                                                                         |
| Emissions test                                                                                                                                                                                        | Compliance | Electromagnetic environment - guidance                                                                                                                                                                                                                                                                                                                                                  |
| RF emissions<br>CISPR 11                                                                                                                                                                              | Group 2    | The <b>RF ABLATOR</b> must emit electromagnetic energy in order to perform its intended function. Nearly electronic equipment may be affected.<br><br>The <b>RF ABLATOR</b> is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| RF emissions<br>CISPR 11                                                                                                                                                                              | Class B    |                                                                                                                                                                                                                                                                                                                                                                                         |
| Harmonic emissions<br>IEC 61000-3-2                                                                                                                                                                   | Class A    |                                                                                                                                                                                                                                                                                                                                                                                         |
| Voltage fluctuations/flicker emissions<br>IEC 61000-3-3                                                                                                                                               | Complies   |                                                                                                                                                                                                                                                                                                                                                                                         |

| Guidance and manufacturer's declaration – electromagnetic immunity                                                                                                                                |                                                                        |                                                                        |                                                                                                                                               |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------|------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|
| The <b>RF ABLATOR</b> is intended for use in the electromagnetic environment specified below. The customer or the user of the <b>RF ABLATOR</b> should assure that is used in such an environment |                                                                        |                                                                        |                                                                                                                                               |
| Immunity test                                                                                                                                                                                     | IEC60601 test level                                                    | Compliance level                                                       | Electromagnetic environment - guidance                                                                                                        |
| Electrostatic discharge (ESD)<br>IEC61000-4-2                                                                                                                                                     | contact $\pm 6$ kV<br>air $\pm 8$ kV                                   | contact $\pm 6$ kV<br>air $\pm 8$ kV                                   | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. |
| Electrical fast transient/burst<br>IEC61000-4-4                                                                                                                                                   | $\pm 2$ kV for power supply lines<br>$\pm 1$ kV for input/output lines | $\pm 2$ kV for power supply lines<br>$\pm 1$ kV for input/output lines | Mains power quality should be that of a typical commercial or hospital environment.                                                           |
| Surge<br>IEC61000-4-5                                                                                                                                                                             | $\pm 1$ kV line(s) to line(s)<br>$\pm 2$ kV line(s) to earth           | $\pm 1$ kV line(s) to line(s)<br>$\pm 2$ kV line(s) to earth           | Mains power quality should be that of a typical commercial or hospital environment.                                                           |

| Guidance and manufacturer's declaration – electromagnetic immunity                                     |                                                                                                                                                                                                                                                                                                                                  |                                                                                                                                                                                                                                                                                                                                  |                                                                                                                                                                                                                                                                                                  |
|--------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Voltage dips, short interruptions and voltage variations on power supply input lines<br>IEC 61000-4-11 | <div>&lt; 5% U<sub>T</sub><br/>(&gt;95% dip in U<sub>T</sub>)<br/>for 0,5 cycle</div> <div>40% U<sub>T</sub><br/>(60% dip in U<sub>T</sub>)<br/>for 5 cycles</div> <div>70% U<sub>T</sub><br/>(30% dip in U<sub>T</sub>)<br/>for 25 cycles</div> <div>&lt; 5% U<sub>T</sub><br/>(&gt;95% dip in U<sub>T</sub>)<br/>for 5 s</div> | <div>&lt; 5% U<sub>T</sub><br/>(&gt;95% dip in U<sub>T</sub>)<br/>for 0,5 cycle</div> <div>40% U<sub>T</sub><br/>(60% dip in U<sub>T</sub>)<br/>for 5 cycles</div> <div>70% U<sub>T</sub><br/>(30% dip in U<sub>T</sub>)<br/>for 25 cycles</div> <div>&lt; 5% U<sub>T</sub><br/>(&gt;95% dip in U<sub>T</sub>)<br/>for 5 s</div> | Mains power quality should be that of a typical commercial or hospital environment.<br>If the user of <b>RF ABLATOR</b> requires continued operation during power mains interruptions, it is recommended that the <b>RF ABLATOR</b> be powered from a uninterruptible power supply or a battery. |
| Power frequency (50/60 Hz) magnetic field<br>IEC 61000-4-8                                             | 3 A/m                                                                                                                                                                                                                                                                                                                            | 3 A/m                                                                                                                                                                                                                                                                                                                            | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment                                                                                                                                                         |

NOTE U<sub>T</sub> is the a.c. mains voltage prior to application of the test level

Guidance and manufacturer's declaration – electromagnetic immunity

The **RF ABLATOR** is intended for use in the electromagnetic environment specified below. The customer or the user of the **RF ABLATOR** should assure that is used in such an environment

| Immunity test                | IEC60601 test level         | Compliance level | Electromagnetic environment - guidance                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|------------------------------|-----------------------------|------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Conducted RF<br>IEC61000-4-6 | 3 Veff<br>150 kHz to 80 MHz | 3 V              | Portable and mobile RF communications equipment should be used no closer to any part of the <b>RF ABLATOR</b> , including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.<br><br><b>Recommended separation distance</b><br><br>$d = 1,2\sqrt{P}$<br>$d = 1,2\sqrt{P}$ from 80 MHz to 800 MHz<br><br>$d = 2,3\sqrt{P}$ from 800 MHz to 2,5 GHz<br><br>where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and $d$ is the recommended separation distance in meters (m).<br><br>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup><br><br>Interference may occur in the vicinity of equipment marked with the following symbol:<br><br> |
| Radiated RF<br>IEC61000-4-3  | 3 V/m<br>80 MHz to 2,5 GHz  | 3 V/m            |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |







- NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.
- NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Le intensità di campo per trasmettitori fissi, come le stazioni base per radiotelefoni (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the **RF ABLATOR** is used exceeds the applicable RF compliance level above, the **RF ABLATOR** should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the **RF ABLATOR**.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

| Recommended separation distances between portable and mobile RF communications equipment and the RF ABLATOR                                                                                                                                                                                                                                                                                                                                                   |                                                                |                                       |                                        |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------|---------------------------------------|----------------------------------------|
| The <b>RF ABLATOR</b> is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the <b>RF ABLATOR</b> can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the <b>RF ABLATOR</b> as recommended below, according to the maximum output power of the communications equipment. |                                                                |                                       |                                        |
| Rated maximum output power of transmitter<br>W                                                                                                                                                                                                                                                                                                                                                                                                                | Separation distance according to frequency of transmitter<br>m |                                       |                                        |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                               | 150 kHz to 80 MHz<br>$d= 1,2\sqrt{P}$                          | 80 MHz to 800 MHz<br>$d= 1,2\sqrt{P}$ | 800 MHz to 2,5 GHz<br>$d= 2,3\sqrt{P}$ |
| 0,01                                                                                                                                                                                                                                                                                                                                                                                                                                                          | 0,12                                                           | 0,12                                  | 0,23                                   |
| 0,1                                                                                                                                                                                                                                                                                                                                                                                                                                                           | 0,38                                                           | 0,38                                  | 0,73                                   |
| 1                                                                                                                                                                                                                                                                                                                                                                                                                                                             | 1,2                                                            | 1,2                                   | 2,3                                    |
| 10                                                                                                                                                                                                                                                                                                                                                                                                                                                            | 3,8                                                            | 3,8                                   | 7,3                                    |
| 100                                                                                                                                                                                                                                                                                                                                                                                                                                                           | 12                                                             | 12                                    | 23                                     |
| For transmitters rated at a maximum output power not listed above, the recommended separation distance $d$ in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.                                                                                                                          |                                                                |                                       |                                        |
| NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.                                                                                                                                                                                                                                                                                                                                                                                             |                                                                |                                       |                                        |
| NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.                                                                                                                                                                                                                                                                                            |                                                                |                                       |                                        |

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